DCM003-12
Ed. 11/2017

ESTRADIOL ELISA

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Ed. 11/2017

DCM003-12

DiaMetra S.r.l. Headquater: Via Calabria 15,
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

DCM003-12
Ed. 11/2017

ESTRADIOL ELISA

for routine analysis

Immunoenzymatic determination of 17 -Estradiol in human serum or plasma.

IVD

LOT

See external label

2°C 8°C

Σ = 96 test

REF DKO003

INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of 17 -Estradiol concentration in human serum or plasma.

Estradiol ELISA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Estradiol (17 -Estradiol) is a sex hormone. It represents the major estrogen in humans. Estradiol has not only a critical impact on reproductive and sexual functioning, but also affects other organs including bone structure. During the reproductive years most Estradiol in women is produced by the ovaries, smaller amounts of Estradiol are also produced by the adrenal cortex. In men, the testes produce Estradiol. In plasma Estradiol is largely bound to sex hormone binding globulin (SHBG), also to albumin, only a fraction is free and biologically active. Serum Estradiol measurement in women reflect primarily the activity of the ovaries. During pregnancy estrogen levels, including Estradiol, rise steadily towards term. Estradiol increases due to placental production. In adult premenopausal women, ovarian production of Estradiol is stimulated by luteinizing hormone (LH) and the follicle-stimulating hormone (FSH) during the menstrual cycle.

In adult women, Estradiol levels are measured in the evaluation of fertility and menstrual irregularities, and to monitor ovarian follicular function during induction of ovulation. In the female, Estradiol acts as a growth hormone for tissue of the reproductive organs. The development of secondary sexual characteristics in women is driven by Estradiol. Estradiol is involved also in man fertility. Estradiol regulates the bone maintenance. Post-menopause women experience an accelerated loss of bone mass due to a relative estrogen deficiency.

Estradiol affects the production of multiple proteins including lipoproteins, binding proteins, and proteins responsible for blood clotting. Estrogens have been found to have neuroprotective function. The Estradiol, for his activities, is involved in some types of cancer such as breast cancer and cancer of the uterine lining. In addition there are several benign gynecologic conditions that are dependent on

estrogen such as endometriosis, leiomyomata uteri, and uterine bleeding.

2. PRINCIPLE

17 -Estradiol (antigen) in the sample competes with the antigenic 17 -Estradiol conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti 17 -Estradiol coated on the microplate (solid phase).

After incubation, the bound/free separation is performed by a simple solid-phase washing. Then the enzyme HRP in the bound-fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H₂SO₄) is added. The colour intensity is inversely proportional to the 17 -Estradiol concentration in the sample.

17 -Estradiol concentration in the sample is calculated through a calibration curve.

3. REAGENT, MATERIAL AND INSTRUMENTATION

3.1. Reagent and material supplied in the kit

1. Calibrators (6 vial)

CAL0 (1 mL)	REF DCE002/0306-0
CAL1 (0.5 mL)	REF DCE002/0307-0
CAL2 (0.5 mL)	REF DCE002/0308-0
CAL3 (0.5 mL)	REF DCE002/0309-0
CAL4 (0.5 mL)	REF DCE002/0310-0
CAL5 (0.5 mL)	REF DCE002/0311-0

2. Control (1 vial, 0.5 mL)

Concentration of Control is indicated on the Certificate of Analysis

REF DCE045/0303-0

3. Conjugate (1 vial, 22 mL)

17 -Estradiol conjugated with horseradish peroxidase (HRP)

REF DCE002/0302-0

4. Coated Microplate (1 microplate breakable)

Anti 17 -Estradiol antibodies adsorbed on microplate

REF DCE002/0303-0

5. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB (0.26 g/L) (avoid any skin contact)

REF DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact)

REF DCE005-0

7. 10X Conc. Wash Solution (1 vial, 50 mL)
Phosphate buffer 0.2M pH 7.4 **REF** DCE054-0

3.2. Reagents necessary not supplied in the kit

Distilled water

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Notes

Store all reagents between 2÷8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of 17 - Estradiol from 20 pg/mL to 2000 pg/mL.
- The clinical significance of 17 -Estradiol determination can be invalidated if the patient was treated with cortisone or natural or synthetic steroids.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.

- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₅)

The Calibrators are ready to use and have the following concentration of 17 -Estradiol:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
pg/mL	0	20	120	300	600	2000

The Calibrators are stable until the expiry date printed on the label. Once opened, the standards are stable six months at 2÷8°C.

6.2. Preparation of the Conjugate

The conjugate is ready to use. Mix gently, for 5 minutes, with a rotating mixer.

Once opened, it is stable six months at 2÷8°C

6.3. Preparation of Wash Solution

Dilute the content of the vial "Conc. Wash Solution 10X" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2÷8°C. In concentrated wash solution is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals, for greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL taking care also to transfer crystals, then mix until crystals are completely dissolved.

6.4. Preparation of the sample

The determination of 17 -Estradiol can be performed in human serum or plasma. Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Avoid repetitive freezing and thawing of samples. Before using, mix gently, for 5 minutes, with a rotating mixer. The Control is ready for use.

6.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C ₀ -C ₅	25 µL		
Sample/Control		25 µL	
Conjugate	200 µL	200 µL	
Incubate 2 h at +37°C. Remove the contents from each well, wash the wells three times with 300 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22÷28°C) for 30 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake gently the microplate. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Estradiol for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Mean Absorbance

Calculate the mean of the absorbance (E_m) for each point of the calibration curve (C₀-C₅) and of each sample

8.2. Calibration curve

Plot the mean value of absorbance (E_m) of the calibrators (C₀-C₅) against concentration. Draw the best-fit curve through the plotted points (es: Four Parameter Logistic).

8.3. Calculation of results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in pg/mL.

9. REFERENCE VALUES

The serum 17 -Estradiol reference values are:

		pg/mL
WOMEN	Follicular phase	30 – 100
	Ovulatory peak	130 – 350
	Luteinic phase	50 – 180
	Menopause	< 60
CHILDREN		< 40
MEN		< 60

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (10x) of three different human sera in one assay. The within assay variability is 9%.

10.1.2. Inter Assay Variation

Between run variation was determined by replicate measurements of three different human sera in different lots. The between assay variability is 10%.

10.2. Accuracy

The dilution test conducted with high concentration samples of 17 β -Estradiol gave an average recovery value (\pm SD) of 95.69% \pm 7.74% with reference to the original concentrations.

The recovery of 120 – 240 – 480 – 960 pg/mL of 17 β -Estradiol added to samples gave an average value (\pm SD) of 101.09% \pm 5.42% with reference to the original concentrations.

10.3. Sensitivity

The lowest detectable concentration of 17 β -Estradiol that can be distinguished from the calibrator zero is 8.68 pg/mL at the 95% confidence limit.

10.4. Specificity

The cross reaction of the antibody calculated at 50% according to Abraham are shown in the table:

Substance	Cross-reactivity
17 β Estradiol	100 %
Estrone	2 %
Estriol	0.39 %
Fulvestrant	0.09 %
Testosterone	0.02 %
Cortisol	< 7x10 ⁻³ %
Progesterone	< 3x10 ⁻⁴ %
DHEA-S	< 1x10 ⁻⁴ %

Important note

Fulvestrant is a chemical compound that is found into the formulation of some drugs used in the treatment of some type of cancers in post-menopausal women; due to its chemical similarity with Estradiol, Fulvestrant molecule can interfere with the assay and lead to an overestimation of Estradiol levels in the sample.

In the case of patients undergoing treatment with Fulvestrant drugs, it is recommended to check the clinical data obtained with Diametra kit with other data for Estradiol quantification in order to verify the interference by Fulvestrant.

10.5. Correlation with a commercial kit

Diametra Estradiol ELISA kit was compared to another commercially available 17 β -Estradiol assay. 16 serum samples were analysed in both test systems.

The linear regression curve is:

$$(17\beta\text{-Estradiol Dia})=1.03*(17\beta\text{-Estradiol Ref})-12.96$$
$$r^2 = 0.996$$

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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