DCM005-12
Ed. 11/2016

DHEA-S ELISA

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Ed. 11/2016

DCM005-12

DiaMetra S.r.l. Headquarter: Via Calabria 15, 20090 SEGRATE (MI) Italy

Tel. +39-02-2139184 - Fax +39-02-2133354

Manufactory: Via Pozzuolo 14, 06038 SPELLO (PG) Italy

Tel. +39-0742-24851 - Fax +39-0742-316197

E-mail: info@diametra.com

DCM005-12
Ed. 11/2016

DHEA-S ELISA

for routine analysis

Direct immunoenzymatic determination of Dehydroepiandrosterone sulfate (DHEA-S) in human serum or plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO005

INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Dehydro-epiandrosterone sulfate (DHEA-S) concentration in human serum or plasma.

DHEA-S ELISA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Dehydroepiandrosterone sulfate (DHEA-S), is a natural steroid hormone found atop of the kidneys in the human body. DHEA-S derived from enzymatic conversion of DHEA in adrenal and extradrenal tissues. DHEA-S is also produced in the gonads, adipose tissue and the brain. It is the most abundant hormone in the human body and it is precursor of all sex steroids.

As most DHEA-S is produced by the zona reticularis of the adrenal, it is argued that there is a role in the immune and stress response. DHEA-S may have more biologic roles. Its production in the brain suggests that is also has a role as a neurosteroid.

Measurement of serum DHEA-S is a useful marker of adrenal androgen synthesis. Abnormally low levels may occur in have been reported in hypoadrenalism, while elevated levels occur in several conditions, e.g. virilizing adrenal adenoma and carcinoma, 21-hydroxylase and 3 -hydroxysteroid dehydrogenase deficiencies and in some cases of female hirsutism. Women with polycystic ovary syndrome tend to have normal or mildly elevated levels of DHEAS. As very little DHEA-S is produced by the gonads, measurement of DHEA-S levels may aid in the localization of androgen source in virilizing conditions.

DHEA-S levels show no diurnal variation.

2. PRINCIPLE OF THE METHOD

The DHEA-S (antigen) in the sample competes with the antigenic DHEA-S conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti DHEA-S coated on the microplate (solid phase). After the incubation, the bound/free separation is performed by a simple solid-phase washing.

Then, the enzyme HRP in the bound-fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H₂SO₄) is added.

The colour intensity is inversely proportional to the DHEA-S concentration of in the sample.

DHEA-S concentration in the sample is calculated through a calibration curve.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 1 mL each, ready to use)

CAL0	REF DCE002/0506-0
CAL1	REF DCE002/0507-0
CAL2	REF DCE002/0508-0
CAL3	REF DCE002/0509-0
CAL4	REF DCE002/0510-0
CAL5	REF DCE002/0511-0

2. Controls (2 vials, 1 mL each)

Control A	REF DCE045/0503A-0
Control B	REF DCE045/0503B-0

Controls Concentration is indicated on the Certificate of Analysis

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5. Coated Microplate (1 breakable microplate)

Antibody anti DHEA-S adsorbed on microplate	REF DCE002/0503-0
---	-------------------

6. TMB Substrate (1 vial, 15 mL)

H ₂ O ₂ -TMB 0.26 g/L (avoid any skin contact)	REF DCE004-0
--	--------------

7. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)	REF DCE005-0
--	--------------

8. 10X Conc. Wash Solution (1 vial, 50 mL)

Phosphate buffer 0.2M pH 7.4	REF DCE054-0
------------------------------	--------------

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Notes

Store all reagents between 20-8°C in the dark.

Open the bag of reagent 5 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit.

Do not remove the adhesive sheets on the unused strips.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants.
- This method allows the determination of Dehydroepiandrosterone Sulphate from 0.1 µg/mL to 10 µg/mL.
- The clinical significance of the determination Dehydroepiandrosterone Sulphate can be invalidated if the patient was treated with cortisone or natural or synthetic steroids.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of

samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of Calibrators (C₀...C₅) and Controls

The Calibrators are ready to use and have the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
µg/mL	0	0.1	0.4	1.0	4.0	10.0

The Calibrators are stable until the expiry date printed on the label. Once opened, the Calibrators are stable for 6 months at 2-8°C.

The Controls are ready to use.

6.2. Preparation of the Sample

The determination of Dehydroepiandrosterone Sulphate can be performed in human plasma as well as in serum of patients.

Store the sample at -20°C if the determination is not performed on the same day of the sample connection. Avoid repetitive freezing and thawing of samples.

Immediately before use, dilute each sample 1:50 with the Serum Diluent (i.e. add 20 µL of sample to 980 µL of Serum Diluent). Mix well.

6.3. Preparation of Wash Solution

Dilute the content of each vial of the "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

In concentrated wash solution is possible to observe the presence of crystals; in this case mix at room temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Diluted Sample/ Control		30 µL	
Calibrators C ₀ -C ₅	30 µL		
Conjugate	100 µL	100 µL	
Incubate at 37°C for 1 hour. Remove the content from each well. Wash the wells 3 times with 300 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22÷28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of DHEA-S for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance

should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Mean Absorbance

Calculate the mean of the absorbances (E_m) for each point of the calibration curve (C₀-C₅) and of each sample.

8.2. Calibration Curve

Plot the values of absorbance (E_m) of the Calibrators (C₀-C₅) against concentration. Draw the best-fit curve through the plotted points (es: Four Parameter Logistic).

8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in µg/mL.

9. REFERENCE VALUES

The serum or plasma Dehydroepiandrosterone Sulphate reference values are:

	WOMAN µg/mL	MAN µg/mL
Newborns	0.9 - 1.8	0.9 - 1.8
Before puberty	0.25 - 1.0	0.25 - 1.0
Adults	0.9 - 3.6	0.9 - 3.6
After menopause	< 0.25 - 1.0	---
Pregnancy	0.25 - 1.8	---

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra-Assay Variation

Within run variation was determined by replicate (16x) the measurements of three different serum samples in one assay. The within assay variability is 7.9%.

10.1.2. Inter-Assay Variation

Between run variation was determined by replicate (20x) the measurement of three different control sera in different lots. The between assay variability is 10.4%.

10.2. Accuracy

The recovery of 0.6 - 1.25 - 2.5 and 5 µg/mL of DHEA-S added to sample gave an average value

(±SD) of 102.87% ± 8.63% with reference to the original concentrations.

The dilution test performed on three sera diluted 2 - 4 and 8 times gave an average value (±SD) of 100.15% ± 9.02%.

10.3. Sensitivity

The lowest detectable concentration of DHEA-S that can be distinguished from the Calibrator 0 is 0.04 µg/mL at the 95% confidence limit.

10.4. Specificity: cross reagent

The cross reaction of the antibody calculated at 50% according to Abraham are shown in the table:

Cross-reagent	Cross reactivity (%)
DHEA-S	100 %
DHEA	100 %
Androstenedione	59 %
Androsterone	16 %
DHT	1,0 %
Testosterone	0,63 %
Estrone	0,3 %
Progesterone	0,27 %
Pregnenolone	0,18 %
17OH Progesterone	0,15 %
11-deoxycortisol	0,08 %
Cortisone	0,013 %
Cortisol	< 0,01 %
17b Estradiol	< 0,01 %
Corticosterone	< 0,01 %
17a Estradiol	< 0,01 %
Cholesterol	< 0,001 %
Estriol	< 0,001 %
Estradiol Sulphate	< 0,001 %
Aldosterone	< 0,001 %
Estradiol-3-Sulphate-17-Glucuronide	< 0,001 %

10.5. Specificity: interfering substances

Interference by Bilirubin, Hemoglobin and Triglycerides has been investigated on Diametra DHEA-S ELISA kit:

Substance	Assayed Conc.	Interference
Bilirubin	0.2 mg/mL	No
Hemoglobin	2 mg/mL	No
Triglycerides	6 mg/mL	No

No interference has been observed with the substances under investigation; following good laboratory practices, it is anyway advised to avoid to use highly lipemic or haemolysed samples.

10.6. Specificity: plasma and SST tube

Interference in plasma and SST (serum separation tube) samples has been evaluated. Serum obtained from the same patient has been used as reference.

Sample	Interference
SST (serum separation tube)	No
EDTA plasma	No
Lithium heparin plasma	No
Sodium heparin plasma	No

No interference has been observed.

10.7. Correlation with RIA

Diametra DHEA-S ELISA kit was compared to another commercially available DHEA-S assay.

42 serum samples have been analysed.

The linear regression curve was calculated:

$$y = 0.93x + 0.28$$

$$r^2 = 0.961$$

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

- Abraham G.E., et al Obstet. Gynecol., 47 (4), 395 (1976)
- Granoff A.B. et al Obstet. Gynecol., 53 (1), 111 (1979)
- Hopper B.R. et al J. Clin. Endocrin. Metab. 40 (3), 458 (1975)
- Winter J.S.D, et al Clin. Obstetric and Gynecol., 21 (1), 67 (1978)
- Cristina Mihaela Ghiciuc C.M et al., Neuroendocrinol Lett 2011; 32(4):475-480

Ed. 11/2016

DCM005-12

DiaMetra S.r.l. Headquarter: Via Calabria 15, 20090 SEGRATE (MI) Italy
Tel. +39-02-2139184
Fax +39-02-2133354

Manufactory: Via Pozzuolo 14, 06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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