DCM006-9
Ed. 01/2015

PROGESTERONE ELISA

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PROGESTERONE ELISA

for routine analysis


Direct immunoenzymatic determination of Progesterone in human serum or plasma

IVD

LOT

See external label

$\Sigma = 96$ tests

REF DKO006

INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Progesterone concentration in human serum or plasma.

Progesterone ELISA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Progesterone is a C-21 steroid hormone involved in the female menstrual cycle, pregnancy (supports *gestation*) and embryogenesis of humans and other species. Progesterone is the major naturally occurring human progestagen.

Progesterone is important for aldosterone (mineralocorticoid) synthesis, as 17-hydroxyprogesterone is for cortisol (glucocorticoid).

Progesterone levels are relatively low in children and postmenopausal women. Adult males have levels similar to those in women during the follicular phase of the menstrual cycle.

In women, progesterone levels are relatively low during the preovulatory phase of the menstrual cycle, rise after ovulation, and are elevated during the luteal phase. If pregnancy occurs, progesterone levels are maintained at luteal levels initially. After delivery of the placenta and during lactation, progesterone levels are very low. The fall in progesterone levels following delivery is one of the triggers for milk production.

Progesterone is produced in the adrenal glands, the gonads (specifically after ovulation in the corpus luteum), the brain, and, during pregnancy, in the placenta.

Progesterone converts the endometrium to its secretory stage to prepare the uterus for implantation. If pregnancy does not occur, progesterone levels will decrease, leading, in the human, to menstruation.

Progesterone belongs to the group of neurosteroids that are found in high concentrations in certain areas in the brain and are synthesized there. Neurosteroids affect synaptic functioning, are neuroprotective, and affect myelination.

Progesterone has multiple effects outside of the reproductive system. Progesterone is thermogenic, it reduces spasm and relaxes smooth muscle. Bronchi are widened and mucus regulated. Progesterone acts as an antiinflammatory agent and regulates the immune response. Progesterone also assists in thyroid function, in bone building by osteoblasts.

Measurement of serum progesterone concentrations have been used in evaluating ovarian function.

2. PRINCIPLE

The Progesterone (antigen) in the sample competes with the antigenic Progesterone conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti Progesterone coated on the microplate (solid phase).

After incubation, the bound/free separation is performed by a simple solid-phase washing.

Then, the enzyme HRP in the bound-fraction reacts with the Substrate (H_2O_2) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H_2SO_4) is added.

The colour intensity is inversely proportional to the Progesterone concentration of in the sample.

Progesterone concentration in the sample is calculated through a calibration curve.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (5 vials, 1 mL each)

CAL0	REF DCE002/0606-0
CAL1	REF DCE002/0607-0
CAL2	REF DCE002/0608-0
CAL3	REF DCE002/0609-0
CAL4	REF DCE002/0610-0

2. Control (1 vial, 1 mL)

Concentration of Control is indicated on the Certificate of Analysis

REF DCE045/0603-0

3. Conjugate (1 vial, 22 mL)

Progesterone conjugated with horseradish peroxidase (HRP)

REF DCE002/0602-0

4. Coated Microplate (1 microplate breakable)

Antibody anti Progesterone coated on the microplate

REF DCE002/0603-0

5. TMB Substrate (1 vial, 15 mL)

H_2O_2 -TMB 0.26 g/L (avoid any skin contact)

REF DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)

REF DCE005-0

7. 10X Conc. Wash Solution (1 vial, 50 mL)
Phosphate buffer 0.2M pH 7.4 **REF** DCE054-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm).

Note

Store all reagents between 20-8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants.
- This method allows the determination of Progesterone from 0.2 ng/mL to 40.0 ng/mL.
- For higher values, for example in pregnancy, dilute the sample; consider the diluting factor when calculating the final result.
- The clinical significance of the Progesterone determination can be invalidated if the patient was treated with cortisone or natural or synthetic steroids.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.

- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrator (C₀...C₄)

Before using leave 5 minutes on a rotary shaker.

The Calibrators are ready to use and have the following concentrations of Progesterone:

	C ₀	C ₁	C ₂	C ₃	C ₄
ng/mL	0	0.2	1.0	8.0	40.0

The Calibrators are stable until the expiry date printed on the label. Once opened, the Calibrators are stable six months at 2-8°C.

6.2. Preparation of the Sample

The determination of Progesterone can be performed in human plasma as well as in serum.

Store reagent at -20°C if the determination is not performed on the same day of the sample connection. Avoid repetitive freezing and thawing of samples.

The Control is ready for use.

6.3. Preparation of the Wash Solution

Dilute the contents of each vial of the "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C. In concentrated wash solution it

is possible to observe the presence of crystals. In this case mix at room temperature until complete dissolution of crystals is observed. For greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer the crystals completely, then mix until the crystals are completely dissolved.

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₄), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Calibrator C ₀ -C ₄	20 µL		
Sample/ Control		20 µL	
Conjugate	200 µL	200 µL	
Incubate at 37°C for 1 hour. Remove the content from each well; wash the wells 3 times with 300 µL of diluted Wash Solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature 22±28°C for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake gently the microplate. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Progesterone for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve (C₀-C₄) and of each sample.

8.2. Calibration curve

Plot the mean value of absorbance (Em) of the Calibrators (C₀-C₄) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in ng/mL.

9. REFERENCE VALUES

The serum or plasma Progesterone reference values are:

		ng/mL	
MEN		< 1.0	
WOMEN	follicular phase	0.1 - 1.4	
	Midluteinic phase	4.0 - 25.0	
	menopause	< 1.0	
	pregnancy	Weeks	
		18-21	53-76
		22-25	60-86
		26-29	71-133
30-33		86-142	
34-37	104-175		
38-41	117-187		

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of

expected values based on the indigenous population where the laboratory works.

Ed. 01/2015

DCM006-9

DiaMetra S.r.l. Headquarter: Via Calabria 15
20090 SEGRATE (MI) Italy
Tel. +39-02-2139184
Fax +39-02-2133354
Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (20x) of two different control sera in one assay. The within assay variability is 4%.

10.1.2. Inter Assay Variation

Between run variation was determined by replicate measurements (10x) of three different control sera in different lots of kit. The between assay variability is 9.3%.

10.2. Accuracy

The recovery of 1.0 - 2.0 - 4.0 - 8.0 ng/mL of Progesterone added to sample gave an average value (\pm SD) of 100.88% \pm 8.29% with reference to the original concentrations.

10.3. Sensitivity

The lowest detectable concentration of Progesterone that can be distinguished from the Calibrator 0 is 0.05 ng/mL at the 95% confidence limit.

10.4. Specificity

The cross reaction of the antibody calculated at 50% according to Abraham are shown in the table:

Progesterone	100 %
Testosterone	0,37 %
17 α OH-progesterone	0,29 %
17 β Estradiolo	0,0013 %
Estrone	0,00053 %
Estriolo	< 0,0001 %
Cortisolo	< 0,0001 %

10.5. Correlation with RIA

Diametra Progesterone ELISA kit was compared to Adaltis Progesterone EIAGen kit commercially available. 31 serum samples were analysed according in both test systems.

The linear regression curve was calculated:

$$(\text{Progest. Diametra}) = 0.97 * (\text{Progest. Adaltis}) + 0.04$$

$$r^2 = 0.887$$

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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