DCM083-5  
Ed. 01/2015

# IAA

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Ed. 01/2015

DCM083-5

**DiaMetra S.r.l. Headquarter:** Via Calabria 15,  
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

**Manufactory:** Via Pozzuolo 14

06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

DCM083-5  
Ed. 01/2015

## IAA

**Quantitative ELISA test for the detection of circulating autoantibodies against human native Insulin in human serum**

for routine analysis

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO083

### INTENDED USE

IAA kit is an enzymatic immunoassay for the quantitative determination of circulating IgG autoantibodies specific to human Insulin.

IAA kit is intended for laboratory use only.

### 1. CLINICAL SIGNIFICANCE

Type 1 diabetes, also known as Insulin-dependent diabetes mellitus (IDDM), results from a chronic autoimmune destruction of the Insulin secreting pancreatic beta cells, probably initiated by exposure of genetically susceptible host to an environmental agent. Autoimmune destruction of beta cells is thought to be completely asymptomatic until 80-90% of the cells are lost. This process may take years to complete and may occur at any time in all ages.

The presence of Insulin autoantibodies (IAA) in patients never treated with Insulin, as opposed to Insulin antibodies (IAb), is an evidence of an ongoing destruction process of pancreatic beta cells in type 1 diabetes.

IAA are particularly important when determining type 1 diabetes risk since their prevalence is significantly elevated in subjects developing the disease in childhood and moreover, they are often the first autoantibodies to be detected before onset of the disease. The prevalence of IAA is inversely correlated with the age of diagnosis.

In type 1 diabetics with recent onset of the disease in the age younger than 5 years IAA can be determined on over 90 % of the patients, whereas in type 1 diabetics over 20 years the prevalence of IAA is greater than 20 %.

The IAA measurement, together with that of antibodies to glutamic acid decarboxylase (GAD65 Ab), protein tyrosine phosphatase-like antigen IA2 and Islet cells antigens (ICA) forms the basis of current strategies for predicting future onset of type 1 diabetes.

In addition to autoantibodies to human Insulin, antibodies induced by Insulin treatment can be frequently found. Such antibodies occur at much higher concentration and are mainly directed against denatured (incorrectly folded) fraction of applied Insulin. This may lead to formation of antigen-antibody-complexes and by this reducing of the activity of the injected Insulin. Measurements are necessary to obtain an optimal disease management when applying Insulin.

### 2. PRINCIPLE

In the first step the antibodies anti Insulin from the diluted sample bind to human recombinant Insulin coated on the microtiter plate. After an incubation of 60 minutes at room temperature (22-28°C), the unbound components are removed by a washing step. In the next step the bound antibodies react with the added anti human IgG antibodies conjugated with horseradish peroxidase (HRP). After 30 minutes at room temperature (22-28°C) the excess of conjugate is removed by another washing step. Then the HRP converts the added colourless TMB Substrate into a blue product. After 15 minutes at room temperature (22-28°C) the enzymatic reaction is stopped by adding the Stop Solution. The absorbance of the resulting yellow product is measured at 450 nm; the obtained optical density is direct proportional to the amount of bound antibodies.

The concentration of antibodies anti Insulin is calculated through a calibration curve.

### 3. REAGENTS, MATERIALS AND INSTRUMENTATION

#### 3.1. Reagent and material supplied in the kit

##### 1. Calibrators (5 vials, 1 mL each)

|      |                   |
|------|-------------------|
| CAL1 | REF DCE002/8307-0 |
| CAL2 | REF DCE002/8308-0 |
| CAL3 | REF DCE002/8309-0 |
| CAL4 | REF DCE002/8310-0 |
| CAL5 | REF DCE002/8311-0 |

##### 2. Controls (2 vials, 1 mL each)

|                  |                   |
|------------------|-------------------|
| Negative Control | REF DCE045/8301-0 |
| Positive Control | REF DCE045/8302-0 |

##### 3. Conjugate (1 vial, 15 mL)

Anti human IgG antibodies conjugated with horseradish peroxidase (HRP)

REF DCE002/8302-0

##### 4. Sample diluent (1 vial, 100 mL)

Phosphate buffer

REF DCE053/8353-0

##### 5. Coated Microplate (1 microplate breakable)

Human recombinant insulin coated on the microplate

REF DCE002/8303-0

##### 6. TMB Substrate (1 vial, 15 mL)

H<sub>2</sub>O<sub>2</sub>-TMB (0.26 g/L) (avoid any skin contact)

REF DCE004/8304-0

7. Stop Solution (1 vial, 15 mL)  
0.25M sulfuric acid (*avoid any skin contact*)  
REF DCE005/8305-0
8. 10X Conc. Wash Solution (1 vial, 100 mL)  
Phosphate buffer REF DCE054/8354-0

### 3.2. Reagents necessary not supplied

Distilled or deionized water.

### 3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

## 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the reagents should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts Proclin 300<sup>R</sup> as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.

## 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.

- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips.** For this purpose, Diametra supplies a separate decontamination reagent for cleaning needles.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of Calibrators (C<sub>1</sub>...C<sub>5</sub>)

The Calibrators are ready for use and have the following concentration:

|      | C <sub>1</sub> | C <sub>2</sub> | C <sub>3</sub> | C <sub>4</sub> | C <sub>5</sub> |
|------|----------------|----------------|----------------|----------------|----------------|
| U/mL | 0.1            | 1              | 5              | 10             | 20             |

The Calibrators are stable until the expiry date printed on the label. Once opened, the Calibrators are stable 6 months at 2-8°C.

### 6.2. Preparation of the Conjugate

The anti human IgG antibodies conjugated with horseradish peroxidase (Conjugate, reactive 3) are ready to use. Once opened the solution is stable up to 4 weeks at 2-8°C.

### 6.3. Preparation of the Wash Solution

Prepare a sufficient amount of washing solution by diluting the 10X Conc. Wash Solution 1:10 with distilled or deionized water. For example, dilute 50 mL of the concentrate Wash Solution with 450 mL of distilled water. The solution should be free of crystals before dilution, otherwise dissolve by warming up to max 37°C. The diluted washing solution can be stored at 2-8°C up to 30 days.

### 6.4. Preparation of the Sample

Blood is taken by venipuncture. After clotting, the serum is separated by centrifugation. Do not use lipaemic or grossly haemolysed serum samples. **Plasma should not be used.** The samples may be kept at 2-8°C up to three days. Long-term storage requires - 20°C. Repeated freezing and thawing should be avoided. For multiple use, aliquot samples and keep them at - 20 °C. **Patients samples have to be diluted 1:101 with sample diluent** (e.g. 5 µL sample + 500 µL sample diluents)

### 6.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.**
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C at maximum 30 days.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>1</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

| Reagent  | Calibrator | Sample/<br>Controls | Blank |
|--|------------|---------------------|-------|
| Diluted sample   |            | 100 µL              |       |
| Calibrator C <sub>1</sub> -C <sub>5</sub>  | 100 µL     |                     |       |
| Controls   |            | 100 µL              |       |
| Cover the plate. Incubate at room temperature (22-28°C) for 1 hour.<br>Remove the contents from each well and wash the wells 3 times with 300 µL of diluted Wash Solution.<br><b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.<br>An insufficient washing will result to poor precision and falsely elevated absorbances. |            |                     |       |

| Conjugate  | 100 µL | 100 µL |        |
|--|--------|--------|--------|
| Cover the plate. Incubate at room temperature (22-28°C) for 30 minutes.<br>Remove the contents from each well and wash the wells 3 times with 300 µL of diluted Wash Solution.<br><b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.<br>An insufficient washing will result to poor precision and falsely elevated absorbances. |        |        |        |
| TMB Substrate  | 100 µL | 100 µL | 100 µL |
| Incubate 15 minutes at room temperature (22-28°C) in the dark  |        |        |        |
| Stop Solution  | 100 µL | 100 µL | 100 µL |
| Shake the microplate gently.<br>Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.  |        |        |        |

## 7. RESULTS

### 7.1. Calibration curve

The calibration curve is established by plotting the mean OD-values of the calibrators 1 - 5 on the ordinate, y-axis, versus their respective IAA-Ab concentrations on the abscissa, x-axis. The IAA Abs concentrations of the controls and the unknown diluted samples are directly read off in U/mL from the measured OD<sub>450</sub> values. There is no further correction for the dilution necessary.

The IAA kit may be used also with Computer Assisted Analysis with software able to use spline smoothing fitting.

Example:

| Sample       | OD (a)<br>450 nm | OD (b)<br>450 nm | OD<br>(mean) | U/mL |
|--------------|------------------|------------------|--------------|------|
| Calibrator 1 | 0.082            | 0.073            | 0.078        | 0.1  |
| Calibrator 2 | 0.288            | 0.226            | 0.257        | 1    |
| Calibrator 3 | 0.824            | 0.720            | 0.772        | 5    |
| Calibrator 4 | 1.804            | 1.700            | 1.752        | 10   |
| Calibrator 5 | 2.697            | 2.607            | 2.652        | 20   |
| Neg. Control | 0.150            | 0.160            | 0.155        | 1.2  |
| Pos. Control | 1.360            | 1.323            | 1.342        | 8.0  |
| Patient 1    | 0.540            | 0.530            | 0.535        | 3.0  |

### Criteria of validation:

Specimens with an OD higher than Calibrator 5 should be diluted further with the sample diluent and the concentration of IAA / IA antibodies should be calculated by applying the dilution factor.

## 7.2. Reference values

| IAA (U/mL) |       |
|------------|-------|
| Negative   | < 2.4 |
| Positive   | ≥ 2.4 |

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

## 8. CHARACTERISTICS

### 8.1. Calibration

The IAA kit is artificially calibrated and concentrations of IAA are therefore expressed in U/mL.

### 8.2. Specificity and Sensitivity

The results available show a diagnostic sensitivity of 77% and a diagnostic specificity of 94% (assayed on 100 sera from healthy subjects). Relative values for sensitivity and specificity refer to corresponding results obtained by the Golden Calibrator IAA RIA. No other possibilities to measure Insulin autoantibodies are known on a routine basis, yet.

### 8.3. Detection Limit

The analytical sensitivity (lower detection limit, blank + 3 SD) is 0.1 U/mL.

The quantification limit (LOQ, blank + 10 SD) is 0.3 U/mL.

### 8.4. Intra and inter-assay variation

| Intra-assay (n=20) |                   |        | Inter-assay (n=5x10) |                   |        |
|--------------------|-------------------|--------|----------------------|-------------------|--------|
| Sample n°          | Mean Conc. (U/mL) | CV (%) | Sample n°            | Mean Conc. (U/mL) | CV (%) |
| 1                  | 4.8               | 7.5    | 4                    | 5.3               | 5.4    |
| 2                  | 9.1               | 2.3    | 5                    | 9.5               | 2.4    |
| 3                  | 18.4              | 4.0    | 6                    | 18.5              | 1.2    |

## 9. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

---

## BIBLIOGRAPHY

1. Hirata Y, H Ishizu, N Ouchi, S Motomura, M Abe, Y Hara,
2. H Wakasugi, I Takahashi, H Sakano, M Tanaka, H Kawao & T Kanasaki: Insulin autoimmunity in a case with spontaneous hypoglycemia; Japan J Diabet 1970, 13: 312-319
3. Palmer JP, CM Asplin, P Clemens, K Lyen, O Tatpati, PK Raghu & TL Paquette: Insulin antibodies in Insulin-dependent diabetes before Insulin treatment; Science 1983, 222:1337-1339
4. Palmer JP, CM Asplin, PK Raghu, P Clemens, K Lyen, O Tatpati, B McKnight, TL Paquette, M Sperling, L Baker & R Guthrie: AntiInsulin antibodies in Insulin-dependent diabetes before Insulin treatment - a new marker for autoimmune beta cell damage?; Pediatr Adolesc Endocrinol 1986, 15:111-116
5. Ziegler, AG, R Ziegler, P Vardi, RA Jackson, JS Soeldner & GS Eisenbarth: Life-table Analysis of Progression to Diabetes of Anti-Insulin Autoantibody-positive Relatives of Individuals with Type 1 Diabetes; Diabetes 1989, 38:1320-1325
6. Williams AJK, PJ Bingley, E Bonifacio, JP Palmer & Eam Gale: A novel Micro-assay for Insulin Autoantibodies; J Autoimmunity 1997; 10:473-478
7. Lindberg B, SA Ivarsson, M Landin-Olsson, G Sundkvist, L Svanberg & A Lernmark: Islet autoantibodies in cord blood from Children who developed Type I (Insulin-dependent) diabetes mellitus before 15 years of age; Diabetologia 1999; 42:181-187
8. Potter KN & T J Wilkins: The molecular specificity of Insulin autoantibodies; Diabetes Metab Res Rev 2000; 16:338-353

Ed. 01/2015

DCM083-5

**DiaMetra S.r.l. Headquarter:** Via Calabria 15,  
20090 SEGRATE (MI) Italy

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**ERROR POSSIBLE CAUSES / SUGGESTIONS****No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

**Too low reaction (too low ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

**Too high reaction (too high ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

**Unexplainable outliers**

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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