DCM087-4
Ed. 01/2015

SHBG

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BIBLIOGRAFIA

- Petra PH. *The plasma sex steroid binding protein (SBP or SHBG). A critical review of recent developments on the structure, molecular biology and function.*J Steroid Biochem Mol Biol. 1991;40(4-6):735-53. Review.
- Selby C. *Sex hormone binding globulin: origin, function and clinical significance.*Ann Clin Biochem. 1990 Nov;27 (Pt 6):532-41. Review.
- Pugeat M, Crave JC, Tourniaire J, Forest MG. *Clinical utility of sex hormone-binding globulin measurement.* Horm Res. 1996;45(3-5):148-55. Review.

Ed. 01/2015

DCM087-4

DiaMetra S.r.l. Headquarter: Via Calabria 15
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

DCM087-4
Ed. 01/2015

SHBG

for routine analysis

Direct immunoenzymatic determination of SHBG in human serum or plasma.

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO087

INTENDED USE

Direct immunoenzymatic colorimetric method for quantitative determination of SHBG concentration in human serum or plasma.

SHBG kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Sex-Hormone-Binding Globulin is a blood transport protein for steroids hormones with an high affinity for testosterone and estradiol.

The SHBG levels are higher in women than in men and increase during pregnancy and in women taking oral contraceptives. Low levels are found in cases of hirsutism, hyperandrogenism and may be encountered in individuals with Cushing's or polycystic ovarian syndrome.

High serum level for men are associated with hypogonadism, hepatic cirrhosis and thyrotoxicosis.

2. PRINCIPLE

SHBG ELISA test is based on the simultaneous binding of human SHBG to two monoclonal antibodies, one immobilized on microwell plates and the other conjugated with horseradish peroxidase (HRP).

After incubation, the bound/free separation is performed by a simple solid-phase washing.

Then the enzyme HRP in the bound-fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H₂SO₄) is added.

The colour intensity is proportional to the SHBG concentration in the sample.

The SHBG concentration in the sample is calculated through a calibration curve..

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (5 vials, 1 mL each)

CAL0	REF DCE002/8706-0
CAL1	REF DCE002/8707-0
CAL2	REF DCE002/8708-0
CAL3	REF DCE002/8709-0
CAL4	REF DCE002/8710-0

2. Control (1 vial, 1 mL, ready to use)

Concentration of Control is indicated on the certificate of Analysis

REF DCE045/8703-0

3. 5X Assay buffer (1 vial, 50 mL)

Hepes buffer; BSA 10 g/L, stabilisers

REF DCE049-0

4. Conjugate (1 vial, 12 mL)

Monoclonal antibody anti SHBG conjugated with horseradish peroxidase (HRP)

REF DCE002/8702-0

5. Coated Microplate (1 breakable microplate)

Monoclonal antibody anti SHBG adsorbed on microplate

REF DCE002/8703-0

6. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB 0.26g/L (avoid any skin contact)

REF DCE004-0

7. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)

REF DCE005-0

8. 50X Conc. Wash Solution (1 vial, 20 mL)

NaCl 45g/L; Tween-20 55g/L

REF DCE006-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm).

Note

The Calibrators contain SHBG in a proteic stabilizing matrix solution.

Store all reagents between 2÷8°C in the dark. Open the bag of reagent 5 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit. Do not remove the adhesive sheets on the strips unutilised.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of SHBG from 2.78 to 277.8 nmol/L.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry dates printed on the labels of the box and of the vials must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment is your responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.

- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated should not be used in the assay. Highly lipemic or haemolysed specimens should similarly not be used
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₄)

The Calibrators have the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄
nmol/L	0	2.78	27.8	111.1	277.8

The Calibrators and Control have already been pre-diluted and are ready to use.

Once opened, the Calibrators are stable 3 months at 2-8°C.

6.2. Preparation of Wash Solution

Dilute the contents of each vial of the buffered Wash Solution Concentrate (50X) with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2÷8°C.

6.3. Preparation of the Sample

Serum or plasma can be used.

Specimen can be stored at 2-8°C for at short time (max two days). For longer storage the specimen should be frozen at -20°C. Avoid repeated freezing and thawing cycles.

Proceed as follows:

1. Dilute the 5X Assay Buffer (1:5):

Dilute 50 mL of assay buffer (reagent 3) with 200 mL of distilled H₂O; for smaller volumes respect the dilution ratio 1:5. Mix gently leaving in a rotating shaker at least 5 minutes.

2. Dilute the samples (1:100):

Dilute 10 µL of sample with 990 µL of diluted assay buffer; for different volumes respect the dilution ratio 1:100.

Mix gently avoiding the use of vortex.

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₄), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C ₀ -C ₅	25 µL		
Diluted Samples/Control		25 µL	
Conjugate	100 µL	100 µL	
Incubate at room temperature (22±28°C) for 1 hours. Remove the content from each well. Wash the wells 3 times with 300 µL of diluted Wash Solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22±28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of SHBG for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Mean Absorbance

Calculate the mean of the absorbancies (E_m) corresponding to the single points of the calibration curve (C₀-C₄) and of each sample.

8.2. Calibration curve

Plot the values of absorbance (E_m) of the Calibrators (C₀-C₄) against concentration. Draw the best-fit curve through the plotted points. (es: Cubic Spline or Four Parameter Logistic).

8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in nmol/L.

9. REFERENCE VALUES

nmol/L	Men	Women
N	102	44
Mean	43	63
Absolute range	15-100	15-120

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Specificity

The antibody specifically recognizes human SHBG. Cross reactivity values have been calculated on a weight/weight basis.

hSHBG (Scripps Laboratories Cat.No.S1724 Lot. 545068)	100 %
hCBG (UCB, Bioproducts Cat. No. i097 Lot. 002/1)	<0.2%
hTBG (Scripps Laboratories Cat.No.T0414 Lot. 780095)	<0.04%

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

- Petra PH. *The plasma sex steroid binding protein (SBP or SHBG). A critical review of recent developments on the structure, molecular biology and function.* J Steroid Biochem Mol Biol. 1991;40(4-6):735-53. Review.
- Selby C. *Sex hormone binding globulin: origin, function and clinical significance.* Ann Clin Biochem. 1990 Nov;27 (Pt 6):532-41. Review.
- Pugeat M, Crave JC, Tourniaire J, Forest MG. *Clinical utility of sex hormone-binding globulin measurement.* Horm Res. 1996;45(3-5):148-55. Review.

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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