DCM019-8
Ed. 03/2017

TOTAL ESTRIOL ELISA

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Ed. 03/2017

DCM019-8

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

DCM019-8
Ed. 03/2017

1

TOTAL ESTRIOLELISA

for routine analysis

Direct immunoenzymatic determination of Total Estriol in human serum or plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO019

INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Total Estriol concentration in human serum or plasma. Total Estriol kit ELISA is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Estriol (also oestriol) is one of the three main estrogens produced by the human body. It is only produced in significant amounts during pregnancy as it is made by the fetus. During pregnancy the production of estriol depends on an intact maternal-placental-fetal unit. Fetal-placental production of estriol leads to a progressive rise in maternal circulating levels reaching a late-gestational peak several orders of magnitude greater than non-pregnant levels. In the maternal circulation, estriol undergoes rapid conjugation in the liver followed by urinary excretion with a half-life of ~20 minutes. Since normal estriol production depends on an intact maternal-placental-fetal circulation and functional fetal metabolism, maternal estriol levels have been used to monitor fetal status during pregnancy, particularly during the third trimester. DHEA is produced by the adrenal cortex of the fetus, this is converted to estriol by the placenta. If levels are abnormally low in a pregnant woman, this may indicate a problem with the development in the child. Levels of estriol in non-pregnant women do not change much after menopause, and levels are not significantly different from levels in men.

2. PRINCIPLE

Total Estriol (antigen) in the sample competes with the antigenic estriol conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti estriol coated on the microplate (solid phase). After incubation, the bound/free separation is performed by a simple solid-phase washing. Then the enzyme HRP in the bound-fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blue color that changes into yellow when the Stop Solution (H₂SO₄) is added. The colour intensity is inversely proportional to the Total Estriol concentration in the sample. Total Estriol concentration in the sample is calculated based on a serie of standards.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

- Calibrators** (5 vials, 1 mL each)

CAL0	REF DCE002/1906-0
CAL1	REF DCE002/1907-0
CAL2	REF DCE002/1908-0
CAL3	REF DCE002/1909-0
CAL4	REF DCE002/1910-0
- Controls** (2 vials, 1 mL each)

Control A	REF DCE045/1903A-0
Control B	REF DCE045/1903B-0

 Concentration of Controls is indicated on the Certificate of Analysis
- Conjugate** (1 vial, 22 mL)
Estriol conjugated with horseradish peroxidase (HRP)
REF DCE002/1902-0
- Coated Microplate** (1 breakable microplate)
Anti estriol antibody adsorbed on microplate
REF DCE002/1903-0
- TMB Substrate** (1 vial, 15 mL)
H₂O₂-TMB 0.26 g/L (avoid any skin contact)
REF DCE004-0
- Stop Solution** (1 vial, 15 mL)
Sulphuric acid 0.15 mol/L (avoid any skin contact)
REF DCE005-0
- 10X Conc. Wash Solution** (1 vial, 50 mL)
Phosphate buffer 0.2M pH 7.4 REF DCE054-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Note

Store all reagents between 2-8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit.

Do not remove the adhesive sheet from the unused strips.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of Total Estriol from 2 ng/mL to 200 ng/mL.
- The clinical significance of Estriol determination can be invalidated if the patient was treated with natural or synthetic steroids.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10

minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators and Controls

The Calibrators are ready to use and have the following concentration of Estriol:

	C ₀	C ₁	C ₂	C ₃	C ₄
ng/mL	0	2	20	80	200

The Controls are ready to use.

Once opened, Calibrators and Controls are stable for 6 months at 2-8°C.

6.2. Preparation of Wash Solution

Dilute the content of each vial of the "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

In concentrated wash solution is possible to observe the presence of crystals; in this case mix at room temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

6.3. Preparation of the Sample

The determination of total Estriol should be performed in human serum or plasma.

Store samples at -20°C if the determination is not performed on the same day of the sample connection. Avoid repetitive freezing and thawing of samples.

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₄), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Sample/ Control		20 µL	
Calibrator C ₀ -C ₄	20 µL		
Conjugate	200 µL	200 µL	
Incubate at 37°C for 1 hour. Remove the content from each well; wash the wells 3 times with 300 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at 22÷28°C for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake gently the microplate. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Total Estriol for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed.

Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established

performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Mean Absorbance

Calculate the mean of the absorbances (E_m) for each point of the calibration curve (C₀-C₄) and of each sample.

8.2. Calibration curve

Plot the values of absorbance (E_m) of the Calibrators (C₀-C₄) against concentration. Draw the best-fit curve through the plotted points(es: Four Parameter Logistic).

8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in ng/mL.

9. REFERENCE VALUES

Serum concentrations of Estriol are included in the following ranges:

Weeks	Median	Range (ng/mL)
17°	18.0	(10 - 27)
18°	25.9	(14 - 51)
19°	39.5	(26 - 52)
20°	40.0	(27 - 53)
21°	45.6	(24 - 66)
22°	39.2	(25 - 58)
23°	56.1	(27 - 70)
24°	56.3	(28 - 75)
25°	64.3	(29 - 84)
26°	68	(41 - 105)
27°	57.4	(41 - 110)
28°	78.0	(38 - 127)
29°	87	(45 - 146)
30°	75	(45 - 160)
31°	88.0	(50 - 170)
32°	90.5	(46 - 175)
33°	100	(60 - 180)
34°	105.6	(60 - 190)
35°	114.2	(65 - 200)
36°	126.0	(74 - 210)
37°	177.0	(90 - 234)
38°	190.0	(101 - 288)
39°	190.0	(102 - 306)
40°	180.0	(60 - 325)
41°	177.5	(95 - 280)

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent

on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay Variation

Within run variation was determined by replicate (16x) the measurement of three different sera in one assay. The within assay variability is 9.9%.

10.1.2. Inter Assay Variation

Between run variation was determined by replicate (10x) the measurement of three different sera in different lots of kit. The between assay variability is 10.3%.

10.2. Accuracy

The recovery of 5.5 - 11 - 22 - 44 ng/mL of Estriol gave an average value (\pm SD) of 103.02% \pm 4.45% with reference to the original concentrations.

The dilution test performed on three sera diluted 2 and 4 times gave an average value (\pm SD) of 107.86% \pm 3.50%

10.3. Sensitivity

The lowest detectable concentration of Total Estriol that can be distinguished from the Calibrator 0 is 1.05 ng/mL at the 95% confidence limit.

10.4. Specificity: cross reagent

The cross reaction of the antibody calculated at 50% according to Abraham are shown in the table:

Cross-reagent	Cross reactivity (%)
Estriol	100 %
16 epi-estriol	10.5 %
15 α OH-estriol	7.0 %
Estriol 3 Sulphate	2.0 %
Estradiol	0.1 %
17 epi-estriol	< 1x10 ⁻² %
Estriol 3 α -Glucuronate	< 1x10 ⁻² %
Estriol 16 α -Glucuronate	< 1x10 ⁻² %
Estrone	< 1x10 ⁻⁴ %

10.5. Specificity: interfering substances

Interference by Bilirubin (conjugated and unconjugated), Hemoglobin and Triglycerides has been investigated on Diametra Total Estriol ELISA kit:

Substance	Assayed Conc.	Interference
Bilirubin (conjugated)	0.2 mg/mL	No
Bilirubin (unconjugated)	0.2 mg/mL	No
Hemoglobin	2 mg/mL	No
Triglycerides	6 mg/mL	No

No interference has been observed with the substances under investigation; following good laboratory practices, it is anyway advised to avoid to use highly lipemic or haemolysed samples.

10.6. Specificity: plasma and SST tube

Interference in plasma and SST (serum separation tube) samples has been evaluated. Serum obtained from the same patient has been used as reference.

Sample	Interference
SST (serum separation tube)	No
EDTA plasma	No
Lithium heparin plasma	No
Sodium heparin plasma	No

No interference has been observed.

10.7. Correlation

The new Diametra Total Estriol ELISA kit was compared to the old Diametra Total Estriol ELISA kit. 35 serum samples was analysed.

The linear regression curve was calculated:

$$Y = 1.02 \cdot X - 1.88$$

$$r^2 = 0.969$$

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

- Fischer-Rasmussen, W., et al Acta Obstet Gynecol. Scand. 60-417 - 420 (1981)
- Truran, P.L., et al Clin. Chem. 28/12, 2393 (1982)
- Vining, R. F., et al J. Clin. Endoc. Metab. 56, 454 (1983)
- Bagger, P.V, et al Acta Obstet Gynecol Scand 60, 187 (1981)
- Osterman, T.M, et al Clin. Chem. 25(5) 716 (1979)
- Wisdom, G.B. Clin. Chem. 22 (8) 1243-1255 (1976)

Ed. 03/2017

DCM019-8

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy
Tel. +39-02-2139184
Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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