DCM056-9
Ed. 01/2015

CA19-9 ELISA

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Ed. 01/2015

DCM056-9

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

DCM056-9
Ed. 01/2015

CA19-9 ELISA

for routine analysis


Direct immunoenzymatic determination of CA19-9 in human serum or plasma.

IVD

LOT

See external label

Σ = 96 tests

REF DKO056

INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of CA19-9 concentration in human serum or plasma.

CA19-9 ELISA is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

CA19-9 represents the most important tumour marker of the mucin type glycoproteins group, that are associated to gastrointestinal cancer pathology.

The immunohistologic distribution of CA19-9 in tissues is consistent with the quantitative determination of higher CA19-9 concentrations in cancer than in normal or inflamed tissues. Recent reports indicate that the serum CA19-9 level is frequently elevated in the serum of subjects with various gastrointestinal malignancies, such as pancreatic, colorectal, gastric and hepatic carcinomas. Together with CEA, elevated CA19-9 is suggestive of gall bladder neoplasm in the setting of inflammatory gall bladder disease. This tumour-associated antigen may also be elevated in some non-malignant conditions.

It has been shown that a persistent elevation in serum CA19-9 value following treatment may be indicative of occult metastatic and/or residual disease. Research studies demonstrate that serum CA19-9 values may have utility in monitoring subjects with the above-mentioned diagnosed malignancies.

A persistently rising serum CA19-9 value may be associated with progressive malignant disease and poor therapeutic response. A declining CA19-9 value may be indicative of a favourable prognosis and good response to treatment.

2. PRINCIPLE

CA19-9 ELISA test is based on simultaneous binding of human CA19-9 to two monoclonal antibodies, one immobilized on microwell plates, the other conjugated with horseradish peroxidase (HRP). After incubation, the separation bound-free is obtained with a simple solid-phase washing. Then, the enzyme in the bound-fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H₂SO₄) is added.

The colour intensity is proportional to CA19-9 concentration in the sample. CA19-9 concentration in the sample is calculated using a calibration curve.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials)

CAL0 (3 mL)	REF DCE002/5606-0
CAL1 (1 mL)	REF DCE002/5607-0
CAL2 (1 mL)	REF DCE002/5608-0
CAL3 (1 mL)	REF DCE002/5609-0
CAL4 (1 mL)	REF DCE002/5610-0
CAL5 (1 mL)	REF DCE002/5611-0

2. Control (1 vial, 1 mL)

The concentration of the Control is lot-specific and is indicated on the Certificate of Analysis

REF DCE045/5603-0

3. Conjugate (1 vial, 12 mL)

Antibody anti CA19-9 conjugated with horseradish peroxidase (HRP)

REF DCE002/5602-0

4. Coated Microplate (1 breakable microplate)

Antibody anti CA19-9 adsorbed on microplate

REF DCE002/5603-0

5. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB 0.26 g/L (avoid any skin contact)

REF DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)

REF DCE005-0

7. 10X Conc. Wash Solution (1 vial, 50 mL)

0.2M Phosphate buffer, pH 7.4

REF DCE054-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Note

Store all reagents at 20±8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until the expiry date of the kit. Do not remove the adhesive sheets on the strips unutilised

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300^R as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of CA19-9 from 0.2 to 240.0 U/mL.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.

- Maximum precision is required for reconstitution and dispensation of reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₅)

The Calibrators are ready to use and have the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
U/mL	0	15	30	60	120	240

Once opened, the Calibrators are stable 6 months at 2-8°C

6.2. Preparation of Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C. In concentrated wash solution is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of the crystals; for greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL taking care also to transfer the crystals, then mix until crystals are completely dissolved.

6.3. Preparation of the Sample

The specimens shall be human serum or plasma; the usual precautions in the collection of venipuncture samples should be observed.

For accurate comparison to established normal values, a fasting morning serum sample should be obtained.

For serum preparation, the blood should be collected in a venipuncture tube without additives or anti-coagulants; then allow the blood to clot and centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, the samples may be stored at -20°C for up to 30 days. Avoid repetitive freezing and thawing.

For sample with concentration over 240 U/mL, dilute the sample with the Calibrator 0.

The Control is ready to use.

6.4. PROCEDURE

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test

results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Sample/ Control		100 µL	
Calibrator C ₀ -C ₅	100 µL		
Incubation 1 h at +37°C. Remove the contents from each well; wash the wells 3 times with 300 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
Conjugate	100 µL	100 µL	
Incubation 1 h at +37°C. Remove the contents from each well; wash the wells 3 times with 300 µL of diluted wash solution. Washing: follow the same indications of the previous point.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature 22÷28°C for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of CA19-9 for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve (C₀-C₅) and of each sample.

8.2. Calibration curve

Plot the mean value of absorbance (Em) of the Calibrators (C₀-C₅) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in U/mL.

9. REFERENCE VALUES

	CA19-9
Healthy people (male and non-pregnant females)	35 U/mL

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay Variation

Within run variation was determined by replicate (10x) the measurement of three different control sera in one assay. The within assay variability is ≤ 3.4%.

10.1.2. Inter Assay Variation

Between run variation was determined by replicate (10x) the measurement of three different control sera with different lots of kit. The between assay variability is ≤ 6.6%.

10.2. Sensitivity

The lowest detectable concentration of CA19-9 that can be distinguished from the Calibrator zero is 0.18 U/mL at the 95 % confidence limit.

10.3. Specificity

The cross-reaction of the CA19-9 assay was assessed by measuring apparent dose of various potentially cross-reagents. The results are shown in the following table:

Antigens	Concentration	% Cross-reactivity
CA19-9	---	100.0%
CA 15-3	1,000 U/mL	3.0%
CA 125	400 U/mL	Not Detected
PSA	100 ng/mL	Not Detected
PAP	1,000 ng/mL	Not Detected
AFP	10,000 ng/mL	Not Detected
CEA	250 ng/mL	Not Detected

Ed. 01/2015

DCM056-9

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

10.4. Dilution

Three patient serum samples containing concentrations of CA19-9 up to 100 IU/mL, were serially diluted with Calibrator zero provided in the kit up to a dilution of 8 and assayed.

The percentage of the mean of recovered dose are 112% (s.d. 3.1%).

10.5. Recovery

Three normal serum samples were spiked with 200, 100, 50, 25 and 0 U/mL of CA19-9 and then assayed.

The percentage of the mean of recovered dose is 92.4% (s.d. 8.8%).

10.6. Correlation

CA19-9 ELISA assay was compared to an established method using Roche Modular; 22 serum samples were assayed with both methods.

A linear regression analysis was performed and the following results were obtained:

(CA19-9 Diametra) = 1.33 (CA19-9 Roche) -18.62

$r^2 = 0.808$

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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