DCM057-6
Ed. 01/2015

tPSA

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BIBLIOGRAFIA

- Christensson A, Laurell CB, Lilja H., *Eur J Biochem*, 194, 755- 63 (1990).
- Watt KW, et. al., *Proc Nat Acad Sci USA*, 83, 3166-70 (1986).
- Chen Z, Prestigiacomo A, Stamey T., *Clin Chem.*, 41, 1273-82 (1995).
- Wild D, *The Immunoassay Handbook.*, Stockton Press (1994) p. 452.
- Junker R, Brandt B, Zechel C, Assmann G, *Clin Chem*, 43, 1588- 94 (1997).
- Prestigiacomo AF, Stamey TA, ' *Physiological variations of serum prostate antigen in the (4-10 ng/ml) range in male volunteers.*' *J. Urol* 1996; 155:1977-80.
- Stamey TA, McNeal JE, Yemoto CM, Sigal BM, Johnstone IM, ' *Biological determinants of cancer progression in men with prostate cancer.*' *JAMA* 1999;281:1395-1400.
- Chen Z, Prestigiacomo A, Stamey T, ' *Purification and characterization of Prostate Specific Antigen (PSA) Complexed to a1 Anticymotrypsin: Potential reference Material for International Standardization of PSA Immunoassays.*' *Clin Chem.* 1995; 41/9:1273-1282.
- Horton GL, Bahnson RR, Datt M, Cphan KM, Catalona WJ and Landenson JH. ' *Differences in values obtained with two assays of Prostate Specific Antigen.*' *J. Urol.* 1988;139:762-72.
- Stenman UH, Leinonen J, Alfthan H, Rannikko S, Tuhkanen K, Alfthan O. ' *A complex between prostate specific antigen and a1-anticymotrypsin is the major form of prostate specific antigen in serum of patients with prostate cancer: assay of complex improves clinical sensitivity for cancer.*' *Cancer Res.* 1991;51:222-26.

Ed. 01/2015

DCM057-6

DiaMetra S.r.l. Headquater: Via Calabria 15,
20090 SEGRATE (MI) Italy
Tel. +39-02-2139184
Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

DCM057-6
Ed. 01/2015

tPSA

Direct immunoenzymatic determination of Total Prostatic Specific Antigen (tPSA) in human serum or plasma

RUO

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO057

INTENDED USE

Diametra tPSA kit is a direct solid phase enzyme immunoassay for quantitative determination of Total Prostatic Specific Antigen (tPSA) in human serum or plasma.

tPSA is intended for laboratory use only.

1. SUMMARY

Prostate Specific Antigen (PSA) is a serine protease with chymotrypsin-like activity (1,2). The protein is a single chain glycoprotein with a molecular weight of 28.4 kDa (3). The name PSA derives from the observation that PSA is a specific antigen of the prostate and it is not found in other tissues.

PSA is found in benign, malignant and metastatic prostate cancer. Since prostate cancer is the second most prevalent form of male malignancy, the detection of elevated PSA levels plays an important role in the early diagnosis. Due to increased sensitivity, serum PSA levels are more useful than prostatic acid phosphatase (PAP) in the diagnosis and management of patients (4).


2. PRINCIPLE OF THE TEST

Immunoenzymometric assay

In this method, PSA calibrators and patient specimen (containing the native PSA antigen) are first added to streptavidin coated wells. Also the biotinylated monoclonal antibody and the enzyme HRP (horseradish peroxidase) labeled antibody are then added (these antibodies are directed against distinct and different epitopes of PSA).

Reaction between the various PSA antibodies and the native PSA of the sample occurs in the microwells without competition or steric hindrance, forming a soluble sandwich complex.

The interaction is illustrated by the following equation:

BtnAb(m) = biotinylated monoclonal antibody (excess quantity)

Ag tPSA = native antigen (variable quantity)

EnzAb = enzyme labeled antibody (excess quantity)

EnzAb -Ag tPSA-BtnAb(m) = antigen-antibodies sandwich complex

k_a = rate constant of association

k_{-a} = rate constant of dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody.

This interaction is illustrated below:

EnzAb -Ag tPSA-BtnAb(m) + StreptavidinC.W.
⇒ **Immobilized complex**

StreptavidinC.W. = Streptavidin immobilized on well
Immobilized complex = Sandwich complex bound to the well

After the incubation period, the sandwich native PSA-antibody-enzyme is separated from unbound conjugate by decantation or aspiration. The enzyme activity in the bound sandwich fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well is quantified by reaction with a suitable substrate to produce color. By utilizing several different calibrators of known PSA values, a dose response curve can be generated from which the unknown PSA concentration in the sample can be ascertained.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

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2. Conjugate (1 vial, 13 mL)

Antibody anti PSA conjugated with horseradish peroxidase (HRP) and Antibody anti PSA biotinylated
REF DCE002/5702-0

3. Coated Microplate (1 breakable microplate)

Microplate coated with streptavidin
REF DCE002/5703-0

4. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB 0.26 g/L (*avoid any skin contact*)
REF DCE004-0

5. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (*avoid any skin contact*)
REF DCE005-0

6. 50X Conc. Wash Solution (1 vial, 20 mL)

NaCl 45 g/L; Tween-20 55 g/L **REF** DCE006-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm).

Notes

Store all reagents at 2-8°C in the dark.

Open the bag of reagent 3 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened it is stable at 2-8°C until the expiry date of the kit.

Do not remove the adhesive sheets on the strips inutilized.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of Proclin 300^R as preservative. Avoid the contact with skin or mucosa.

- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of tPSA from 5.0 ng/mL to 100 ng/mL.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators

The Calibrators are ready to use, are calibrated against the 1st IS 96/670, and have the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
ng/mL	0	5	10	25	50	100

Once opened, the Calibrators are stable 6 months at 2-8°C.

6.2. Preparation of Wash Solution

Dilute the contents of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

6.3. Preparation of the Sample

The specimens can be human serum or plasma; the usual precautions in the collection of venipuncture samples should be observed.

For accurate comparison to established normal values, a fasting morning serum sample should be obtained.

To obtain the serum, the blood should be collected in a venipuncture tube without additives or anti-coagulants; allow the blood to clot; centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, they may be stored at -20°C for up to 30 days. Avoid repetitive freezing and thawing.

When assayed in duplicate, 0,050 mL of the specimen is required.

Patient specimens with tPSA concentrations above 100 ng/mL may be diluted (for example 1:10 or higher) with normal female serum (tPSA < 5 ng/mL) and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor.

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagents	Calibrator	Campione	Blank
Calibrator C ₀ -C ₅	25 µL		
Sample		25 µL	
Conjugate	100 µL	100 µL	
Incubate at room temperature (22-28°C) for 30 minutes. Remove the content from each well; wash the wells 3 times with 300 µL of diluted Wash Solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22-28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at levels in the low, medium and high ranges of the dose response curve for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Note

The optical densities (OD) of some calibrators and samples may be higher than 2.0, and, in such a case, they could be out of the measurement range of the microplate reader. It is therefore necessary, for OD higher than 2.0, to perform a reading at 405 nm (= wavelength of peak shoulder) in addition to 450 nm (peak wavelength) and 620 (reference filter for the subtraction of interferences due to the plastic). For microplate readers unable to read the plate at 3 wavelengths at the same time, it is advisable to proceed as follows:

- read the microplate at 450 nm and at 620 nm.
- read again the plate at 405 nm and 620 nm.

- find out the wells whose OD at 450 nm are higher than 2.0
- select the corresponding OD read at 405 nm and multiply these values at 405 nm by the conversion factor 3.0 (where $OD_{450nm} / OD_{405nm} = 3.0$), that is: $OD_{450nm} = OD_{405nm} \times 3.0$.

Warning: The conversion factor 3.0 is suggested only. For better accuracy, the user is advised to calculate the conversion factor specific for his own reader.

8.2. Data Processing - Automated Method

Use the 4 parameters logistic – preferred – or the smoothed cubic spline function as calculation algorithm.

NOTE: If computer controlled data reduction is used to calculate the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

8.3. Data Processing - Manual Method

A dose response curve is used to ascertain the concentration of tPSA in unknown specimens.

- Record the OD obtained from the printout of the microplate reader as outlined in Example 1.
- Draw on a linear graph paper a dose response curve (DRC) using the mean OD for each duplicate calibrator versus the corresponding tPSA concentration in ng/mL.
- Draw the best-fit curve through the plotted points.
- To determine the concentration of tPSA for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/mL) in the horizontal axis of the graph.

EXAMPLE 1				
Sample ID	Well number	OD	Mean OD	Value ng/mL
CAL 0	A1	0.019	0.019	0
	B1	0.019		
CAL 1	C1	0.279	0.276	5
	D1	0.273		
CAL 2	E1	0.563	0.561	10
	F1	0.559		
CAL 3	G1	1.248	1.213	25
	H1	1.179		
CAL 4	A2	2.051	1.999	50
	B2	1.947		
CAL 5	C2	2.892	2.833	100
	D2	2.775		
patient	E2	1.186	1.142	23.6
	F2	1.089		

The data presented in Example 1 are for illustration only and should not be used in lieu of a dose response curve prepared with each assay.

9. REFERENCE VALUES

Healthy males are expected to have values below 4 ng/mL (4).

It is important to keep in mind that establishment of a values range which can be expected to found by a given method for a population of "normal"-persons is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually. It is recommended that each laboratory establishes its own normal and pathological ranges of seric tPSA.

10. INTERPRETATION OF THE ASSAY

PSA is elevated in benign prostrate hypertrophy (BPH). **Clinically an elevated PSA value alone is not of diagnostic value as a specific test for cancer** and should only be used in conjunction with other clinical manifestations and diagnostic procedures (prostate biopsy). Free PSA determinations may be helpful in regard to the discrimination of BPH and prostrate cancer conditions (5).

11. PERFORMANCE AND CHARACTERISTICS

11.1. Precision

The within and between assay precision of the tPSA kit were determined by analysis on three different levels of control sera. The number, mean value, standard deviation (SD) and coefficient of variation for each of these control sera are presented in Table 1 and Table 2.

TABLE 1				
Within Assay Precision (Values in ng/mL)				
Sample	N	Mean	SD	C.V.
Level 1	20	0.7	0.05	7.1%
Level 2	20	4.5	0.20	4.4%
Level 3	20	28.3	1.07	3.7%

TABLE 2				
Between Assay Precision (Values in ng/mL)				
Sample	N	Mean	SD	C.V.
Level 1	10	0.8	0.09	11.3%
Level 2	10	4.3	0.25	5.8%
Level 3	10	27.5	1.42	5.2%

11.2. Correlation

The tPSA kit was compared to a reference Elisa method. Biological specimens with normal and elevated concentrations were assayed. The total number of such specimens was 241. The least square regression equation and the correlation

coefficient were computed for the tPSA kit in comparison with the reference method. The data obtained is displayed in Table 3.

Method	Mean	Least Square Regression Analysis	Correlation Coefficient
EIAgen tPSA kit	5.62	$y = -0.0598 + 0.98(x)$	0.987
Reference method	5.57		

11.3. Sensitivity

The tPSA kit has a sensitivity of 0.012 ng. This is equivalent to a sample containing 0.5 ng/mL tPSA concentration.

11.4. Specificity

No interference was detected with the performance of tPSA Elisa upon addition of massive amounts of the following substances to a human serum pool.

Substance	Concentration
Acetylsalicylic Acid	100 µg/mL
Ascorbic Acid	100 µg/mL
Caffeine	100 µg/mL
CEA	10 µg/mL
AFP	10 µg/mL
CA-125	10000 U/mL
hCG	1000 IU/mL
hLH	10 IU/mL
hTSH	100 mIU/mL
hPRL	100 µg/mL

12. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

- Christensson A, Laurell CB, Lilja H., *Eur J Biochem*, 194, 755- 63 (1990).
- Watt KW, et. al., *Proc Nat Acad Sci USA*, 83, 3166-70 (1986).
- Chen Z, Prestigiacomo A, Stamey T., *Clin Chem*, 41, 1273-82 (1995).
- Wild D, *The Immunoassay Handbook*, Stockton Press (1994) p. 452.
- Junker R, Brandt B, Zechel C, Assmann G, *Clin Chem*, 43, 1588- 94 (1997).
- Prestigiacomo AF, Stamey TA, ' *Physiological variations of serum prostate antigen in the (4-10 ng/ml) range in male volunteers.*' *J. Urol* 1996; 155:1977-80.
- Stamey TA, McNeal JE, Yemoto CM, Sigal BM, Johnstone IM, ' *Biological determinants of cancer progression in men with prostate cancer.*' *JAMA* 1999;281:1395-1400.
- Chen Z, Prestigiacomo A, Stamey T, ' *Purification and characterization of Prostate Specific Antigen (PSA) Complexed to a1 Anticymotrypsin: Potential reference Material for International Standardization of PSA Immunoassays.*' *Clin Chem*. 1995; 41/9:1273-1282.
- Horton GL, Bahnson RR, Datt M, Cphan KM, Catalona WJ and Landenson JH. ' *Differences in values obtained*

with two assays of Prostate Specific Antigen.' *J. Urol*. 1988;139:762-72.

- Stenman UH, Leinonen J, Alfthan H, Rannikko S, Tuhkanen Kand Alfthan O. ' *A complex between prostate specific antigen and a1-anticymotrypsin is the major form of prostate specific antigen in serum of patients with prostate cancer: assay of complex improves clinical sensitivity for cancer.*' *Cancer Res*. 1991;51:222-26.

Ed. 01/2015

DCM057-6

DiaMetra S.r.l. Headquarter: Via Calabria 15,
20090 SEGRATE (MI) Italy
Tel. +39-02-2139184
Fax +39-02-2133354

Manufactory: Via Pozzuolo 14,
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

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- Christensson A, Laurell CB, Lilja H., *Eur J Biochem*, 194, 755- 63 (1990).
- Watt KW, et. al., *Proc Nat Acad Sci USA*, 83, 3166-70 (1986).
- Chen Z, Prestigiacomo A, Stamey T., *Clin Chem.*, 41, 1273-82 (1995).
- Wild D, *The Immunoassay Handbook.*, Stockton Press (1994) p. 452.
- Junker R, Brandt B, Zechel C, Assmann G, *Clin Chem*, 43, 1588- 94 (1997).
- Prestigiacomo AF, Stamey TA, ' *Physiological variations of serum prostate antigen in the (4-10 ng/mL) range in male volunteers.*' *J. Urol* 1996;155:1977-80.
- Stamey TA, McNeal JE, Yemoto CM, Sigal BM, Johnstone IM, ' *Biological determinants of cancer progression in men with prostate cancer.*' *JAMA* 1999;281:1395-1400.
- Chen Z, Prestigiacomo A, Stamey T, ' *Purification and characterization of Prostate Specific Antigen (PSA) Complexed to a1-Anticymotrypsin: Potential reference Material for International Standardization of PSA Immunoassays.*' *Clin Chem.* 1995;41/9:1273-1282.
- Horton GL, Bahnson RR, Datt M, Cfhan KM, Catalona WJ and Landenson JH. ' *Differences in values obtained with two assays of Prostate Specific Antigen.*' *J. Urol.* 1988;139:762-72.
- Stenman UH, Leinonen J, Alfthan H, Rannikko S, Tuhkanen Kand Alfthan O. ' *A complex between prostate specific antigen anda1-anticymotrypsin is the major form of prostate specific antigenin serum of patients with prostate cancer: assay of*

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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