

兔抗 CCND3(Phospho-Thr283) 多克隆抗体

中文名称：兔抗 CCND3(Phospho-Thr283) 多克隆抗体

英文名称： Anti-CCND3(Phospho-Thr283) rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

抗 原： CCND3(Phospho-Thr283)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene belongs to the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event.
Applications:	WB
Name of antibody:	CCND3(Phospho-Thr283)
Immunogen:	Peptide sequence around phosphorylation site of Threonine 283 (T-S-T(p)-P-T) derived from Human Cyclin D3.

Full name:	cyclin D3
Synonyms :	0
SwissProt:	P30281
WB Predicted band size:	33 kDa
WB Positive control:	Mouse heart tissue lysates
WB Recommended dilution:	500-1000

