

兔抗 DNAJC14 多克隆抗体

中文名称：兔抗 DNAJC14 多克隆抗体

英文名称：Anti-DNAJC14 rabbit polyclonal antibody

别名：DnaJ heat shock protein family (Hsp40) member C14; DNAJ; HDJ3; LIP6; DRIP78

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DNAJC14

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Regulates the export of target proteins, such as DRD 1, from the endoplasmic reticulum to the cell surface .
Applications:	ELISA, IHC
Name of antibody:	DNAJC14
Immunogen:	Fusion protein of human DNAJC14
Full name:	DnaJ heat shock protein family (Hsp40) member C14
Synonyms:	DNAJ; HDJ3; LIP6; DRIP78
SwissProt:	Q6Y2X3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil

IHC Recommend dilution:

50-300

