

兔抗 CCNG1 多克隆抗体

中文名称：兔抗 CCNG1 多克隆抗体

英文名称： Anti-CCNG1 rabbit polyclonal antibody

别名： CCNG

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CCNG1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The eukaryotic cell cycle is governed by cyclin-dependent protein kinases (CDKs) whose activities are regulated by cyclins and CDK inhibitors. The protein encoded by this gene is a member of the cyclin family and contains the cyclin box. The encoded protein lacks the protein destabilizing (PEST) sequence that is present in other family members. Transcriptional activation of this gene can be induced by tumor protein p53. Two transcript variants encoding the same protein have been identified for this gene.
Applications:	ELISA, IHC

Name of antibody:	CCNG1
Immunogen:	Synthetic peptide of human CCNG1
Full name:	cyclin G1
Synonyms :	CCNG
SwissProt:	P51959
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and human tonsil
IHC Recommend dilution:	25-100



