

兔抗 DNASE1L1 多克隆抗体

中文名称：兔抗 DNASE1L1 多克隆抗体

英文名称： Anti-DNASE1L1 rabbit polyclonal antibody

别名： deoxyribonuclease 1 like 1; XIB; G4.8; DNL1L; DNASEX; DNAS1L1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DNASE1L1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a deoxyribonuclease protein that shows high sequence similarity to DNase I. The encoded protein is localized to the endoplasmic reticulum and modified by N-linked glycosylation. Alternate transcriptional splice variants encoding the same protein have been observed.
Applications:	ELISA, WB, IHC
Name of antibody:	DNASE1L1
Immunogen:	Fusion protein of human DNASE1L1
Full name:	deoxyribonuclease 1 like 1
Synonyms:	XIB; G4.8; DNL1L; DNASEX; DNAS1L1

SwissProt:	P49184
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	34 kDa
WB Positive control:	Human heart tissue lysate
WB Recommended dilution:	200-1000



