

兔抗 DNASE1L2 多克隆抗体

中文名称： 兔抗 DNASE1L2 多克隆抗体

英文名称： Anti-DNASE1L2 rabbit polyclonal antibody

别名： deoxyribonuclease I-like 2; DNASE1L2

相关类别： 一抗

抗原： DNASE1L

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Divalent cation-dependent acid DNA endonuclease involved in the breakdown of the nucleus during cornification of epidermal keratinocytes. May play an immune role by eliminating harmful DNA released into the extracellular environment by damaged epidermal cells.
Applications:	ELISA, IHC
Name of antibody:	DNASE1L2
Immunogen:	Synthetic peptide of human DNASE1L2
Full name:	deoxyribonuclease I-like 2
Synonyms:	DNASE1L2
SwissProt:	Q92874

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human colorectal cancer
IHC Recommend dilution:	25-100



