

兔抗 CCNI 多克隆抗体

中文名称：兔抗 CCNI 多克隆抗体

英文名称： Anti-CCNI rabbit polyclonal antibody

别名： cyclin I; CYI; CYC1; CCNI1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CCNI

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The protein encoded by this gene belongs to the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. This cyclin shows the highest similarity with cyclin G. The transcript of this gene was found to be expressed constantly during cell cycle progression.

Applications:	ELISA, IHC
Name of antibody:	CCNI
Immunogen:	Fusion protein of human CCNI
Full name:	cyclin I
Synonyms:	CYI; CYC1; CCNI1
SwissProt:	Q14094
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human colorectal cancer
IHC Recommend dilution:	50-200



