

兔抗 DNM1L 多克隆抗体

中文名称：兔抗 DNM1L 多克隆抗体

英文名称：Anti-DNM1L rabbit polyclonal antibody

别名：dynamin 1 like; DLP1; DRP1; DVLP; EMPF; OPA5; EMPF1; DYMPLE; HDYNIV

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DNM1L

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the dynamin super family of GTPases. The encoded protein mediates mitochondrial and peroxisomal division, and is involved in developmentally regulated apoptosis and programmed necrosis. Dysfunction of this gene is implicated in several neurological disorders, including Alzheimer's disease. Mutations in this gene are associated with the autosomal dominant disorder, encephalopathy, lethal, due to defective mitochondrial and peroxisomal fission (EMPF). Alternative splicing results in multiple transcript variants encoding different isoforms.

	orms. [provided by RefSeq, Jun 2013]
Applications:	ELISA, WB, IHC
Name of antibody:	DNM1L
Immunogen:	Fusion protein of human DNM1L
Full name:	dynamamin 1 like
Synonyms:	DLP1; DRP1; DVLP; EMPF; OPA5; EMPF1; DYMPLE; H DYNIV
SwissProt:	O00429
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human colorectal cancer
IHC Recommend dilution:	200-300
WB Predicted band size:	82 kDa
WB Positive control:	Raji cell and Mouse brain tissue lysates
WB Recommended dilution:	500-2000



