

## 兔抗 DNM3 多克隆抗体

中文名称：兔抗 DNM3 多克隆抗体

英文名称： Anti-DNM3 rabbit polyclonal antibody

别名： Dyna III

抗原： DNM3

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of a family of guanosine triphosphate (GTP)-binding proteins that associate with microtubules and are involved in vesicular transport. The encoded protein functions in the development of megakaryocytes. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	DNM3
<b>Immunogen:</b>	Fusion protein of human DNM3
<b>Full name:</b>	dynamin 3
<b>Synonyms :</b>	Dyna III

<b>SwissProt:</b>	Q9UQ16
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	98 kDa
<b>WB Positive control:</b>	Mouse brain tissue
<b>WB Recommended dilution:</b>	500-2000

