

兔抗 CD1D 多克隆抗体

中文名称：兔抗 CD1D 多克隆抗体

英文名称： Anti-CD1D rabbit polyclonal antibody

别名： R3; CD1A

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CD1D

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a divergent member of the CD1 family of transmembrane glycoproteins, which are structurally related to the major histocompatibility complex (MHC) proteins and form heterodimers with beta-2-microglobulin. The CD1 proteins mediate the presentation of primarily lipid and glycolipid antigens of self or microbial origin to T cells. The human genome contains five CD1 family genes organized in a cluster on chromosome 1.
Applications:	ELISA, WB, IHC
Name of antibody:	CD1D

Immunogen:	Fusion protein of human CD1D
Full name:	CD1d molecule
Synonyms :	R3; CD1A
SwissProt:	P15813
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100
WB Predicted band size:	38 kDa
WB Positive control:	Jurkat cell
WB Recommended dilution:	500-2000



