

兔抗 DOK4 多克隆抗体

中文名称： 兔抗 DOK4 多克隆抗体

英文名称： Anti-DOK4 rabbit polyclonal antibody

别名： docking protein 4 IRS5; IRS-5

储存： 冷冻 (-20℃)

抗原： DOK4

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	DOK proteins are enzymatically inert adaptor or scaffolding proteins. They provide a docking platform for the assembly of multimolecular signaling complexes. DOK4 functions in RET-mediated neurite outgrowth and plays a positive role in activation of the MAP kinase pathway (By similarity). Putative link with downstream effectors of RET in neuronal differentiation. May be involved in the regulation of the immune response induced by T-cells.
Applications:	ELISA, IHC
Name of antibody:	DOK4
Immunogen:	Full length fusion protein

Full name:	docking protein 4
Synonyms:	IRS5; IRS-5
SwissProt:	Q8TEW6
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human esophagus cancer; Human prostate cancer



