

兔抗 DOP1A 多克隆抗体

中文名称：兔抗 DOP1A 多克隆抗体

英文名称： Anti-DOP1A rabbit polyclonal antibody

别名： DOP1 leucine zipper like protein A; DOP1; DOPEY1;

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	May be involved in protein traffic between late Golgi and early endosomes.
Applications:	ELISA, IHC
Name of antibody:	DOP1A
Immunogen:	Synthetic peptide of human DOP1A
Full name:	DOP1 leucine zipper like protein A
Synonyms:	DOP1; DOPEY1; KIAA1117; dJ202D23.2
SwissProt:	Q5JWR5
IHC positive control:	Human liver cancer and Human ovarian cancer
IHC Recommend dilution:	25-100

