

兔抗 DPF2 多克隆抗体

中文名称：兔抗 DPF2 多克隆抗体

英文名称： Anti-DPF2 rabbit polyclonal antibody

别名： REQ; UBID4; ubi-d4

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： DPF2

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the d4 domain family, characterized by a zinc finger-like structural motif. This protein functions as a transcription factor which is necessary for the apoptotic response following deprivation of survival factors. It likely serves a regulatory role in rapid hematopoietic cell growth and turnover. This gene is considered a candidate gene for multiple endocrine neoplasia type I, an inherited cancer syndrome involving multiple parathyroid, enteropancreatic, and pituitary tumors.
Applications:	ELISA, WB

Name of antibody:	DPF2
Immunogen:	Fusion protein of human DPF2
Full name:	D4, zinc and double PHD fingers family 2
Synonyms :	REQ; UBID4; ubi-d4
SwissProt:	Q92785
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	44 kDa
WB Positive control:	K562 and Jurkat cells
WB Recommended dilution:	200-1000

