

兔抗 DPPA2 多克隆抗体

中文名称：兔抗 DPPA2 多克隆抗体

英文名称： Anti-DPPA2 rabbit polyclonal antibody

别名： CT100; PESCRG1; ECAT15-2

抗原： DPPA2

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This protein binds to target gene promoters, including NKX2-5 and SYCE1, but not GATA4, and may be involved in the maintenance of the active epigenetic status of these genes. Expressed in embryonic stem cells. No expression is seen in 5 months embryo, mesenchymal stem cells, embryonic fibrocytes and adult tissues. It contains 1 SAP domain.
Applications:	ELISA, IHC
Name of antibody:	DPPA2
Immunogen:	Fusion protein of human DPPA2
Full name:	developmental pluripotency associated 2

Synonyms :	CT100; PESCRG1; ECAT15-2
SwissProt:	Q7Z7J5
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	25-100

