

兔抗 CD34 多克隆抗体

中文名称：兔抗 CD34 多克隆抗体

英文名称：Anti-CD34 rabbit polyclonal antibody

储存：冷冻（-20℃）

抗原：CD34

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene may play a role in the attachment of stem cells to the bone marrow extracellular matrix or to stromal cells. This single-pass membrane protein is highly glycosylated and phosphorylated by protein kinase C. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	CD34
Immunogen:	Synthetic peptide of human CD34
Full name:	CD34 molecule
SwissProt:	P28906
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and human liver cancer

IHC Recommend dilution:	10-50
-------------------------	-------



