

兔抗 DPT 多克隆抗体

中文名称：兔抗 DPT 多克隆抗体

英文名称： Anti-DPT rabbit polyclonal antibody

别名： TRAMP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DPT

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Dermatopontin is an extracellular matrix protein with possible functions in cell-matrix interactions and matrix assembly. The protein is found in various tissues and many of its tyrosine residues are sulphated. Dermatopontin is postulated to modify the behavior of TGF-beta through interaction with decorin.
Applications:	ELISA, IHC
Name of antibody:	DPT
Immunogen:	Fusion protein of human DPT
Full name:	dermatopontin
Synonyms:	TRAMP

SwissProt:	Q07507
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	20-100



