

ADA2 抗原（重组蛋白）

中文名称：ADA2 抗原（重组蛋白）

英文名称：ADA2 Antigen (Recombinant Protein)

别名：adenosine deaminase 2; PAN; ADGF; CECR1; IDGFL; SNEDS; VAIHS

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 312-511 amino acids of human ADA2

技术规格

Full name:	adenosine deaminase 2
Synonyms:	PAN; ADGF; CECR1; IDGFL; SNEDS; VAIHS
Swissprot:	Q9NZK5
Gene Accession:	BC051755
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of a subfamily of the adenosine deaminase protein family. The encoded protein is one of two ade

nosine deaminases found in humans, which regulate levels of the signaling molecule, adenosine. The encoded protein is secreted from monocytes undergoing differentiation and may regulate cell proliferation and differentiation. This gene may be responsible for some of the phenotypic features associated with cat eye syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]