

## 兔抗 DPYSL2(Ab-509) 多克隆抗体

中文名称：兔抗 DPYSL2(Ab-509) 多克隆抗体

英文名称：Anti-DPYSL2(Ab-509) rabbit polyclonal antibody

别名：DRP2; N2A3; CRMP2; DRP-2; ULIP2; CRMP-2; DHPRP2; ULIP-2

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：DPYSL2(Ab-509)

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Plays a role in neuronal development and polarity, as well as in axon growth and guidance, neuronal growth cone collapse and cell migration. Necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton. May play a role in endocytosis.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	DPYSL2(Ab-509)
<b>Immunogen:</b>	Synthesized non-phosphopeptide derived from human CRMP-2 around the phosphorylation site of threonine 509 (S-V-T(p)-P-K).
<b>Full name:</b>	dihydropyrimidinase-like 2

<b>Synonyms :</b>	DRP2; N2A3; CRMP2; DRP-2; ULIP2; CRMP-2; DHPRP2; ULI P-2
<b>SwissProt:</b>	Q16555
<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	62 kDa
<b>WB Positive control:</b>	HT-29 cells treated with heat shock lysate
<b>WB Recommended dilution:</b>	500-3000



