

兔抗 DPYSL2(Phospho-Ser522) 多克隆抗体

中文名称：兔抗 DPYSL2(Phospho-Ser522) 多克隆抗体

英文名称： Anti-DPYSL2(Phospho-Ser522) rabbit polyclonal antibody

别名： DRP2; N2A3; CRMP2; DRP-2; ULIP2; CRMP-2; DHPRP2; ULIP-2

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： DPYSL2(Phospho-Ser522)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the collapsin response mediator protein family. Collapsin response mediator proteins form homo- and hetero-tetramers and facilitate neuron guidance, growth and polarity. The encoded protein promotes microtubule assembly and is required for Sema3A-mediated growth cone collapse, and also plays a role in synaptic signaling through interactions with calcium channels.
Applications:	WB
Name of antibody:	DPYSL2(Phospho-Ser522)
Immunogen:	Peptide sequence around phosphorylation site of Serine 522 (K-T-S(p)-P-A) derived from Human CRMP-2.

Full name:	dihydropyrimidinase-like 2
Synonyms :	DRP2; N2A3; CRMP2; DRP-2; ULIP2; CRMP-2; DHPRP2; ULIP-2
SwissProt:	Q16555
WB Predicted band size:	62 kDa
WB Positive control:	Huvec cells lysates
WB Recommended dilution:	500-1000

