

兔抗 DPYSL4 多克隆抗体

中文名称：兔抗 DPYSL4 多克隆抗体

英文名称： Anti-DPYSL4 rabbit polyclonal antibody

别名： CRMP3, DRP-4, ULIP4

相关类别： 一抗

抗原： DPYSL4

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Collapsin response mediator protein 3 is an enzyme that in humans is encoded by the CRMP3 gene. It is necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton. And plays a role in axon guidance, neuronal growth cone collapse and cell migration. Homotetramer, and heterotetramer with CRMP1, DPYSL2, DPYSL3 or DPYSL5. Interacts with PLEXA1. Phosphorylated up on DNA damage, probably by ATM or ATR.
Applications:	ELISA, WB, IHC
Name of antibody:	DPYSL4

Immunogen:	Synthetic peptide of human DPYSL4
Full name:	dihydropyrimidinase-like 4
Synonyms :	CRMP3, DRP-4, ULIP4
SwissProt:	O14531
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human cervical cancer and Human ovarian cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	62 kDa
WB Positive control:	Human fetal brain tissue
WB Recommended dilution:	500-1000



