

兔抗 DPYSL5 多克隆抗体

中文名称: 兔抗 DPYSL5 多克隆抗体

英文名称: Anti-DPYSL5 rabbit polyclonal antibody

别 名: CRAM; CRMP5; Ulip6; CRMP-5

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: DPYSL5

反应种属: Human, Mouse, Rat

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the CRMP (collapsi ng response mediator protein) family thought to be involved in neural development. Antibodies to the e ncoded protein were found in some patients with n eurologic symptoms who had paraneoplastic syndro me. A pseudogene of this gene is found on chromo some 11. Multiple alternatively spliced variants, enco ding the same protein, have been identified.
Applications:	ELISA, WB
Name of antibody:	DPYSL5
Immunogen:	Synthetic peptide of human DPYSL5



Full name:	dihydropyrimidinase-like 5
Synonyms:	CRAM; CRMP5; Ulip6; CRMP-5
SwissProt:	Q9BPU6
ELISA Recommended dilution:	500-5000
WB Predicted band size:	62 kDa
WB Positive control:	Human fetal brain tissue
WB Recommended dilution:	200-1000

