

兔抗 CD80 多克隆抗体

中文名称：兔抗 CD80 多克隆抗体

英文名称：Anti-CD80 rabbit polyclonal antibody

别名：B7; BB1; B7-1; B7.1; LAB7; CD28LG; CD28LG1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CD80

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a membrane receptor that is activated by the binding of CD28 or CTLA-4. The activated protein induces T-cell proliferation and cytokine production. This protein can act as a receptor for a denovirus subgroup B and may play a role in lupus neuropathy.
Applications:	ELISA, IHC
Name of antibody:	CD80
Immunogen:	Synthetic peptide of human CD80
Full name:	CD80 molecule
Synonyms :	B7; BB1; B7-1; B7.1; LAB7; CD28LG; CD28LG1

SwissProt:	P33681
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human gastric cancer
IHC Recommend dilution:	25-100

