

## 兔抗 DRGX 多克隆抗体

中文名称：兔抗 DRGX 多克隆抗体

英文名称：Anti-DRGX rabbit polyclonal antibody

别名：dorsal root ganglia homeobox; DRG11; PRRXL1

相关类别：一抗

抗原：DRGX

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Transcription factor required for the formation of correct projections from nociceptive sensory neurons to the dorsal horn of the spinal cord and normal perception of pain.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DRGX
<b>Immunogen:</b>	Synthetic peptide of human DRGX
<b>Full name:</b>	dorsal root ganglia homeobox
<b>Synonyms:</b>	DRG11; PRRXL1
<b>SwissProt:</b>	A6NNA5
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human breast cancer
<b>IHC Recommend dilution:</b>	50-200

