

兔抗 DSC1 多克隆抗体

- 中文名称: 兔抗 DSC1 多克隆抗体
- 英文名称: Anti-DSC1 rabbit polyclonal antibody
- 别名: CDHF1; DG2/DG3
- 抗 原: DSC1
- 储存: 冷冻(-20℃)
- 宿 主: Rabbit
- 相关类别: 一抗
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a calcium-depende nt glycoprotein that is a member of the desmocollin su bfamily of the cadherin superfamily. These desmosomal family members, along with the desmogleins, are found primarily in epithelial cells where they constitute the ad hesive proteins of the desmosome cell-cell junction and are required for cell adhesion and desmosome formatio n. The desmosomal family members are arranged in tw o clusters on chromosome 18, occupying less than 650 kb combined. Alternative splicing results in two transcri
	pt variants encoding distinct isoforms.



Applications:	ELISA, IHC
Name of antibody:	DSC1
Immunogen:	Synthetic peptide of human DSC1
Full name:	desmocollin 1
Synonyms:	CDHF1; DG2/DG3
SwissProt:	Q08554
IHC positive control:	Human gastric cancer and human brain
IHC Recommend dilution:	50-200

