

兔抗 CD151 多克隆抗体

中文名称：兔抗 CD151 多克隆抗体

英文名称： Anti-CD151 rabbit polyclonal antibody

别名： GP27, MER2, RAPH, SFA1, PETA-3, TSPAN24

抗原： CD151

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell surface glycoprotein that is known to complex with integrins and other transmembrane 4 superfamily proteins. It is involved in cellular processes including cell adhesion and may reg

	ulate integrin trafficking and/or function. This protein enhances cell motility, invasion and metastasis of cancer cells. Multiple alternatively spliced transcript variants that encode the same protein have been described for this gene.
Applications:	ELISA, IHC
Name of antibody:	CD151
Immunogen:	Fusion protein of human CD151
Full name:	Cluster of Differentiation 151 antigen
Synonyms :	GP27, MER2, RAPH, SFA1, PETA-3, TSPAN24
SwissProt:	P48509
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200

