

兔抗 CD200R1 多克隆抗体

中文名称：兔抗 CD200R1 多克隆抗体

英文名称： Anti-CD200R1 rabbit polyclonal antibody

别名： OX2R; MOX2R; CD200R; HCRTR2

相关类别： 一抗

储存： 冷冻 (-20℃)

抗原： CD200R1

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a receptor for the OX-2 membrane glycoprotein. Both the receptor and substrate are cell surface glycoproteins containing two immunoglobulin-like domains. This receptor is restricted to the surfaces of myeloid lineage cells and the receptor-substrate interaction may function as a myeloid downregulatory signal. Mouse studies of a related gene suggest that this interaction may control myeloid function in a tissue-specific manner. Alternative splicing of this gene results in multiple transcript variants.
Applications:	ELISA, IHC

Name of antibody:	CD200R1
Immunogen:	Fusion protein of human CD200R1
Full name:	CD200 receptor 1
Synonyms :	OX2R; MOX2R; CD200R; HCRTR2
SwissProt:	Q8TD46
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human prostate cancer
IHC Recommend dilution:	50-200

