

兔抗 DSG1 多克隆抗体

中文名称：兔抗 DSG1 多克隆抗体

英文名称：Anti-DSG1 rabbit polyclonal antibody

别名：DG1; DSG; CDHF4; PPKS1; SPPK1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DSG1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

Desmosomes are cell-cell junctions between epithelial, myocardial and certain other cell types. Desmoglein 1 is a calcium-binding transmembrane glycoprotein component of desmosomes in vertebrate epithelial cells. Currently, three desmoglein subfamily members have been identified and all are members of the cadherin cell adhesion molecule superfamily. These desmoglein gene family members are located in a cluster on chromosome 18. The protein encoded by this gene has been identified as the autoantigen of the autoimmune skin blistering disease pemphigus foliaceus.

Applications:	ELISA, WB
Name of antibody:	DSG1
Immunogen:	Synthetic peptide of human DSG1
Full name:	desmoglein 1
Synonyms :	DG1; DSG; CDHF4; PPKS1; SPPK1
SwissProt:	Q02413
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	114 kDa
WB Positive control:	Human tongue cancer tissue
WB Recommended dilution:	200-1000

