

兔抗 DSG2 多克隆抗体

中文名称：兔抗 DSG2 多克隆抗体

英文名称： Anti-DSG2 rabbit polyclonal antibody

别名： HDGC; CDHF5; ARVC10; ARVD10; CMD1BB

相关类别： 一抗

储存： 冷冻 (-20℃)

抗原： DSG2

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Desmosomes are cell-cell junctions between epithelial, myocardial, and certain other cell types. This gene product is a calcium-binding transmembrane glycoprotein component of desmosomes in vertebrate epithelial cells. Currently, three desmoglein subfamily members have been identified and all are members of the cadherin cell adhesion molecule superfamily. These desmoglein gene family members are located in a cluster on chromosome 18. This second family member is expressed in colon, colon carcinoma, and other simple and stratified epithelial-derived cell lines.

Applications:	ELISA, IHC
Name of antibody:	DSG2
Immunogen:	Fusion protein of human DSG2
Full name:	desmoglein 2
Synonyms :	HDGC; CDHF5; ARVC10; ARVD10; CMD1BB
SwissProt:	Q14126
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and human gastric cancer
IHC Recommend dilution:	50-200

