

## 兔抗 CD226 多克隆抗体

中文名称：兔抗 CD226 多克隆抗体

英文名称：Anti-CD226 rabbit polyclonal antibody

别名：PTA1; DNAM1; DNAM-1; TLISA1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CD226

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a glycoprotein expressed on the surface of NK cells, platelets, monocytes and a subset of T cells. It is a member of the Ig-superfamily containing 2 Ig-like domains of the V-set. The protein mediates cellular adhesion of platelets and megakaryocytic cells to vascular endothelial cells. The protein also plays a role in megakaryocytic cell maturation.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CD226
<b>Immunogen:</b>	Synthetic peptide of human CD226
<b>Full name:</b>	CD226 molecule

<b>Synonyms :</b>	PTA1; DNAM1; DNAM-1; TLISA1
<b>SwissProt:</b>	Q15762
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	50-200

