

兔抗 DSTN 多克隆抗体

中文名称：兔抗 DSTN 多克隆抗体

英文名称： Anti-DSTN rabbit polyclonal antibody

别名： ADF; ACTDP; HEL32; bA462D18.2

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： DSTN

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The product of this gene belongs to the actin-binding proteins ADF family. This family of proteins is responsible for enhancing the turnover rate of actin in vivo. This gene encodes the actin depolymerizing protein that severs actin filaments (F-actin) and binds to actin monomers (G-actin). Two transcript variants encoding distinct isoforms have been identified for this gene.
Applications:	ELISA, IHC
Name of antibody:	DSTN
Immunogen:	Fusion protein of human DSTN
Full name:	destrin (actin depolymerizing factor)

Synonyms :	ADF; ACTDP; HEL32; bA462D18.2
SwissProt:	P60981
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human liver cancer
IHC Recommend dilution:	50-200

