

HCCS 抗原（重组蛋白）

中文名称：HCCS 抗原（重组蛋白）

英文名称：HCCS Antigen (Recombinant Protein)

别名：MLS; CCHL; MCOPS7; LSDMCA1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-268 amino acids of human HCCS

技术规格：

Full name:	holocytochrome c synthase
Synonyms:	MLS; CCHL; MCOPS7; LSDMCA1
Swissprot:	P53701
Gene Accession:	BC001691
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is an enzyme that covalently links a heme group to the apoprotein of cytochrome c. Defects in this gene are a cause of microphthalmia syndromic type 7 (MCOPS7). Three transcript variants encoding the same protein have been found for this gene.