

兔抗 FLT1(Phospho-Tyr1048) 多克隆抗体

中文名称：兔抗 FLT1(Phospho-Tyr1048) 多克隆抗体

英文名称： Anti-FLT1(Phospho-Tyr1048) rabbit polyclonal antibody

别名： FLT; FLT-1; VEGFR1; VEGFR-1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： FLT1(Phospho-Tyr1048)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a member of the vascular endothelial growth factor receptor (VEGFR) family. VEGFR family members are receptor tyrosine kinases (RTKs) which contain an extracellular ligand-binding region with seven immunoglobulin (Ig)-like domains, a transmembrane segment, and a tyrosine kinase (TK) domain within the cytoplasmic domain. This protein binds to VEGFR-A, VEGFR-B and placental growth factor and plays an important role in angiogenesis and vasculogenesis. Expression of this receptor is found in vascular endothelial cells, placental trophoblast cells and peripheral blood monocytes. Multiple tra

	nscript variants encoding different isoforms have been found for this gene. Isoforms include a full-length transmembrane receptor isoform and shortened, soluble isoforms. The soluble isoforms are associated with the onset of pre-eclampsia.[provided by RefSeq, May 2009]
Applications:	IHC
Name of antibody:	FLT1(Phospho-Tyr1048)
Immunogen:	Peptide sequence around phosphorylation site of tyrosine 1048 (D-I-Y(p)-K-N) derived from Human VEGFR1.
Full name:	fms-related tyrosine kinase 1
Synonyms :	FLT; FLT-1; VEGFR1; VEGFR-1
SwissProt:	P17948
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100

