

兔抗 FLT1 多克隆抗体

中文名称：兔抗 FLT1 多克隆抗体

英文名称：Anti-FLT1 rabbit polyclonal antibody

别名：FLT; FLT-1; VEGFR1; VEGFR-1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：FLT1

反应种属：Human Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Receptor for VEGF, VEGFB and PGF. Has a tyrosine-protein kinase activity. The VEGF-kinase ligand/receptor signaling system plays a key role in vascular development and regulation of vascular permeability. Isoform SFlt1 may have an inhibitory role in angiogenesis.
Applications:	WB
Name of antibody:	FLT1
Immunogen:	Synthetic peptide of human FLT1
Full name:	fms-related tyrosine kinase 1
Synonyms :	FLT; FLT-1; VEGFR1; VEGFR-1

SwissProt:	P17948
WB Predicted band size:	151 kDa
WB Positive control:	Molt4, Hela, MCF-7 and PC12 cells
WB Recommended dilution:	500-1000

