

兔抗 FLT3LG 多克隆抗体

中文名称：兔抗 FLT3LG 多克隆抗体

英文名称： Anti-FLT3LG rabbit polyclonal antibody

别名： FL; FLT3L

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： FLT3LG

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Dendritic cells (DCs) provide the key link between innate and adaptive immunity by recognizing pathogens and priming pathogen-specific immune responses. FLT3LG controls the development of DCs and is particularly important for plasmacytoid DCs and CD8 (see MIM 186910)-positive classical DCs and their CD103 (ITGAE; MIM 604682)-positive tissue counterparts (summary by Sathaliyawala et al., 2010 [PubMed 20933441]).
Applications:	WB
Name of antibody:	FLT3LG
Immunogen:	Synthesized peptide derived from internal of human Flt3 I

	igand.
Full name:	fms-related tyrosine kinase 3 ligand
Synonyms :	FL; FLT3L
SwissProt:	P49771
WB Predicted band size:	26 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

