

## 兔抗 FLT4 多克隆抗体

中文名称: 兔抗 FLT4 多克隆抗体

英文名称: Anti-FLT4 rabbit polyclonal antibody

级别: PCL; FLT-4; FLT41; LMPH1A; VEGFR3; VEGFR-3

抗原: FLT4

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a tyrosine kinase receptor for vascular endothelial growth factors C and D. The protein is thought to be involved in lymphangiogenesis and maintenance of the lymphatic endothelium. Mutations in this gene cause hereditary lymphedema type IA.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FLT4
<b>Immunogen:</b>	Synthetic peptide of human FLT4
<b>Full name:</b>	fms related tyrosine kinase 4
<b>Synonyms:</b>	PCL; FLT-4; FLT41; LMPH1A; VEGFR3; VEGFR-3
<b>SwissProt:</b>	P35916

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	30-150



