

兔抗 FMN2 多克隆抗体

中文名称：兔抗 FMN2 多克隆抗体

英文名称：Anti-FMN2 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

宿 主：Rabbit

抗 原：FMN2

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Formin homology (FH) domain proteins (see FMN1; MIM 136535) play a role in cytoskeletal organization and /or establishment of cell polarity.
Applications:	ELISA, IHC
Name of antibody:	FMN2
Immunogen:	Synthetic peptide of human FMN2
Full name:	formin 2
SwissProt:	Q9NZ56
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain and Human colon cancer
IHC Recommend dilution:	25-100

