

## 兔抗 FMNL1 多克隆抗体

- 中文名称：兔抗 FMNL1 多克隆抗体
- 英文名称：Anti-FMNL1 rabbit polyclonal antibody
- 别名：FMNL; FHOD4; KW-13; C17orf1; C17orf1B
- 相关类别：一抗
- 储存：冷冻 (-20℃)
- 抗原：FMNL1
- 宿主：Rabbit
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a formin-related protein. Formin-related proteins have been implicated in morphogenesis, cytokinesis, and cell polarity. An alternative splice variant has been described but its full length sequence has not been determined.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	FMNL1
<b>Immunogen:</b>	Fusion protein of human FMNL1
<b>Full name:</b>	formin-like 1
<b>Synonyms :</b>	FMNL; FHOD4; KW-13; C17orf1; C17orf1B
<b>SwissProt:</b>	O95466

<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	20-100
<b>WB Predicted band size:</b>	122 kDa
<b>WB Positive control:</b>	A172 cells
<b>WB Recommended dilution:</b>	500-2000



