

IQCD 抗原（重组蛋白）

中文名称：IQCD 抗原（重组蛋白）

英文名称：IQCD Antigen (Recombinant Protein)

别名：DRC10; CFAP84; 4933433C09Rik

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-200 amino acids of human IQCD

技术规格：

Full name:	IQ motif containing D
Synonyms:	DRC10; CFAP84; 4933433C09Rik
Swissprot:	Q96DY2
Gene Accession:	BC013151
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	IQCD contains 1 IQ domain. The function of the IQCD protein is not known.