

## KANSL1L 抗原（重组蛋白）

中文名称：KANSL1L 抗原（重组蛋白）

英文名称： KANSL1L Antigen (Recombinant Protein)

别名： MSL1v2; C2orf67

储存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 477-676 amino acids of human KANSL1L

技术规格：

<b>Full name:</b>	KAT8 regulatory NSL complex subunit 1 like
<b>Synonyms:</b>	MSL1v2; C2orf67
<b>Swissprot:</b>	A0AUZ9
<b>Gene Accession:</b>	BC032837
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The function of this protein remains unknown.