

兔抗 FOSL1 多克隆抗体

中文名称： 兔抗 FOSL1 多克隆抗体

英文名称： Anti-FOSL1 rabbit polyclonal antibody

别名： FRA; FRA1; fra-1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： FOSL1

反应种属： Human, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation.
Applications:	ELISA, WB
Name of antibody:	FOSL1
Immunogen:	Synthetic peptide of human FOSL1
Full name:	FOS-like antigen 1

Synonyms :	FRA; FRA1; fra-1
SwissProt:	P15407
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	29 kDa
WB Positive control:	231 cells
WB Recommended dilution:	500-2000

