

兔抗 FOSL2 多克隆抗体

中文名称：兔抗 FOSL2 多克隆抗体

英文名称： Anti-FOSL2 rabbit polyclonal antibody

别名： FOS like 2, AP-1 transcription factor subunit; FRA2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： FOSL2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation.
Applications:	ELISA, WB
Name of antibody:	FOSL2
Immunogen:	Fusion protein of human FOSL2
Full name:	FOS like 2, AP-1 transcription factor subunit
Synonyms:	FRA2

SwissProt:	P15408
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	35 kDa
WB Positive control:	Mouse adrenal gland tissue lysate
WB Recommended dilution:	500-2000

