

兔抗 FOXA1 多克隆抗体

中文名称：兔抗 FOXA1 多克隆抗体

英文名称：Anti-FOXA1 rabbit polyclonal antibody

别名：HNF3A; TCF3A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FOXA1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver.
Applications:	ELISA, WB, IHC
Name of antibody:	FOXA1
Immunogen:	Synthetic peptide of human FOXA1
Full name:	forkhead box A1

Synonyms :	HNF3A; TCF3A
SwissProt:	P55317
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	49 kDa
WB Positive control:	Human prostate tissue and MCF-7 cells
WB Recommended dilution:	200-1000



