

## 兔抗 FOXB2 多克隆抗体

中文名称：兔抗 FOXB2 多克隆抗体

英文名称：Anti-FOXB2 rabbit polyclonal antibody

别名：bA159H20.4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FOXB2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

FOX (Forkhead box) proteins are a family of transcription factors that play important roles in regulating the expression of genes involved in cell growth, proliferation, differentiation, and longevity. Many FOX proteins are important to embryonic development. The defining feature of FOX proteins is the forkhead box, a sequence of 80 to 100 amino acids forming a motif that binds to DNA. This forkhead motif is also known as the winged helix due to the butterfly-like appearance of the loops in the protein structure of the domain. Forkhead genes are a subgroup of the helix-turn-helix class of proteins. The FOXF2 gene encodes fo

	rkhead box F2, one of many human homologues of the Drosophila melanogaster transcription factor forkhead. FOXF2 is expressed in lung and placenta.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FOXB2
<b>Immunogen:</b>	Synthetic peptide of human FOXB2
<b>Full name:</b>	Forkhead box B2
<b>Synonyms :</b>	bA159H20.4
<b>SwissProt:</b>	Q5VYV0
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	10-50

